



Datasheet
V1.42.20

Contents

1. System Overview	5
1.1 Architecture	5
1.2 Applications	6
1.3 Imaging features	6
1.4 Standard features	7
1.5 Language support	8
2. Ergonomics.....	8
2.1 Keyboard	9
2.2 Touch screen.....	9
2.3 Image display screen	9
2.4 Wheels.....	9
2.5 Touch gestures	9
2.6 System boot-up.....	10
2.7 Comments	10
2.8 Bodymark.....	10
2.9 Peripherals.....	10
2.10 Dimensions and Weight	10
2.11 Electrical Power.....	10
2.12 Operating Environment	10
2.13 Storage & Transportation Environment	10
3. Transducers	11
3.1 Transducer Technology	11
3.2 Transducer types.....	11
3.3 Transducer selection	11
4. Advanced Imaging controls.....	18
4.1 VFusion	18
4.2 VSpeckle	18
4.3 VTissue	18
4.4 3D/4D.....	18
4.4.1 3D/4D HQ Grad.....	18
4.4.2 Inversion mode	18
4.4.3 Magic Cut(optional)	19
4.4.4 Niche view.....	19
4.4.5 Free View(optional).....	19
4.4.6 3D Smart Face(optional).....	19
4.4.7 VNavIn (optional).....	19
4.5 Tissue Doppler (TD).....	19
4.6 Tissue Velocity Imaging (TVI).....	20
4.7 Stress Echo(optional).....	20
4.8 Strain Imaging(optional).....	20
4.9 Tissue Velocity M mode(TVM) (optional)	20
4.10 Multi-angle M mode.....	20

4.11 Curved M mode	21
4.12 VAid (optional)	21
5. Imaging modes	21
5.1 2D Imaging	21
5.2 Harmonic Imaging	22
5.3 M mode	22
5.4 Color Doppler mode	22
5.5 Power Doppler mode	23
5.6 Pulsed Wave (PW) Doppler	23
5.7 Continuous Wave Doppler (CW)	24
5.8 Elastography imaging	24
5.9 3D/4D	24
5.10 PView (optional)	25
5.11 TView	26
5.12 Auto optimization	26
5.13 Live Track (including angle free flow)(optional)	26
5.14 Tissue Doppler Imaging (TD)	26
5.15 Tissue Velocity Imaging (TVI)	26
5.16 Tissue Velocity M mode (TVM) (optional)	27
6. Touch Panel Interface	27
6.1 2D mode	27
6.2 M Mode	28
6.3 CF mode	28
6.4 PW/CW mode	28
6.5 3D mode	29
6.6 4D mode	29
7. System Feature	30
7.1 Display modes	30
7.2 Display annotation	30
7.3 Simple User Operation Interface	31
7.4 Cineloop	31
7.5 Quick save feature	31
7.6 ECG (optional)	31
7.7 Archive	31
7.8 Report	32
7.9 Connectivity	32
7.10 Probes/application	32
7.11 Safety Conformance	33
8. Measurement and Analysis	33
8.1 Measurement in different modes	33
8.1.1 Generic Measurement in 2D mode	33
8.1.2 Generic Measurement in CM mode	34
8.1.3 Generic Measurement in M mode	34
8.1.4 Measurement in PW mode	34
8.2 Measurement in different applications	35

8.2.1 Abdominal Measurement.....	35
8.2.2 Small Part Measurement.....	35
8.2.3 Vessel Measurement	35
8.2.4 Gynecology Measurement.....	35
8.2.5 Urology Measurement	35
8.2.6 Pediatric Measurement	35
8.2.7 Obstetrics Measurement	35
8.2.8 Cardiac Measurement	35
8.2.9 Auto NT (Nuchal Translucency) measurement(optional).....	36
8.2.10 Auto IMT (Intima-Media Thickness) measurement.....	36
8.2.11 Live IMT (Intima-Media Thickness) measurement (optional).....	36
8.2.12 Auto IT (Intracranial translucency) measurement(optional)	36
8.2.13 Auto Follicle(2D/3D)(3D optional)	36
8.2.14 Smart 3D Volume Measurement(optional).....	36
8.2.15 VAim OB measurement(optional).....	36
8.2.16 VAim Hip measurement(optional)	37
8.2.17 VAim Follicle(2D)measurement (optional)	37
8.2.18 VAid Breast(optional).....	37
8.2.19 VAid Thyroid(optional).....	37

Ultrasound System Specifications

Extremely portable and exceptional performance ES-R500 meets all your clinical needs by:

- Unmatched image quality
- All ranges of features, functions, and probes
- Flexible and customized simple workflow
- Powered by Artificial intelligent technologies
- Extremely flexible and compact design with 23.8 inch monitor

1. System Overview

1.1 Architecture

- ES-R500 brings a confident diagnostic experience with the extraordinary processing power of our breakthrough VLucid⁺ Platform, to deliver superior image quality, thanks to its exceptional intelligent architecture
- The new generation VLucid⁺ platform is capable of processing multiple data streams simultaneously
- The new 12 bit, low noise, digital circuitry, with up to 280db dynamic range has improved 2D performance



and increased Doppler sensitivity

- Directional-enhanced information compiling for more tissue detail and reduction of angle-generated artifacts
- Next generation adaptive image processing for noise and artifact reduction that improves tissue presentation and edge definition
- Fully independent, triplex multiple mode operation for easy in Doppler procedures
- Multi-processors allow simultaneous mode changes and support for advanced system functionality
 - VLuminous Flow provides the color Doppler flow innovatively in a 3D view with excellent sensitivity, which can help understand the structure of blood flow and small vessels intuitively
 - Sync ROI enables the width of 2D scan area is synchronized with the CF ROI, which effectively improves the frame rate
 - Zscore analysis, provide a new way for fetal heart evaluation
 - Support to export 3D data for 3D printer(optional)

- Support multiple DICOM server configuration(optional)
- Background transfer, supports background export without interrupting the actual scan
- Foot switch(optional)
- VReport, a customer-centric tool for report templates design, makes the whole report procedure more smooth and individual(optional)
- Customized user interface, allows user to change the position of buttons on the touch screen, also the size of 'probe&app' UI window is adjustable
- VWork, an intelligent feature, which enables users to configure workflows for every application scenario. This leads to easy and effective adherence to a department protocol and saves operation time to a great extent

1.2 Applications

- Abdomen
- Obstetric
- Gynecology
- Cardiology
- Urology
- Vascular
- TCD
- Small Parts
- Pediatrics

1.3 Imaging features

- 2D grayscale imaging
- Harmonic
- VFusion
- VSpeckle
- VTissue
- Auto imaging optimization
- Live Track(optional)
- Easy Comparative Function to compare previous exam
- M mode
- Color Doppler
- Power Doppler
- Pulse wave Doppler
- Multi Doppler(optional)
- PWV(optional)
- Simultaneous 2D and M mode
- Duplex 2D/PW Doppler
- Triplex 2D/Color/PW Doppler
- High PRF pulsed wave Doppler
- Continuous wave Doppler
- Zoom
- FULL screen imaging to enlarge imaging size
- Dual real time imaging without compromising imaging size
- PView for panoramic imaging (optional)
- TView for trapezoidal imaging
- Elastography imaging
- Needle Enhancement
- SGC (Scanline gain compensation)
- Cardiac Quantification(optional)

- 2D auto follicle
- Free 3D (optional)
- 3D/4D imaging
- HQ 3D/4D
- HQ Silhouette (optional)
- Tomographic display (MCUT)
- Inversion mode
- Magic Cut(optional)
- Free View(optional)
- Niche view
- Light Lab(optional)
- Color 3D(optional)
- ECG(optional)
- Tissue Doppler (TD) mode
- Tissue Velocity Imaging (TVI)
- Tissue Velocity M (TVM) (optional)
- Stress echo(optional)
- Strain imaging(optional)
- Auto NT(Nuchal translucency) (optional)
- Auto IT, automatic measurement of Intracranial translucency(optional)
- VLuminous flow, a feature which shows the blood flow in a 3-D view with excellent sensitivity
- VFlow,adaptive color flow filter to increase the sensitivity of blood flow
- VAid Breast(optional)
- VAid Thyroid(optional)
- VAid HRI(optional)
- Color M-mode
- Curved M mode, user can draw any curved sample line freely and get corresponding results

- Multi-line Angular M-Mode, Up to 4 sample lines
- Sync B/C width, the width of B mode interest area is always be the same with the CF mode
- Live IMT, display intima-media thickness in real time(optional)
- VAim(optional)
(Vino for OB, Follicle, Hip, Pelvic, LEVA)
- 3D Smart Face, an intelligent tool for fetal face optimization(optional)
- VNavIn, a tool that navigates inside the 3D volume data and projects an inside-out perspective image that displays the inner most structures like virtual endoscopy(optional)

1.4 Standard features

- Up to 1500 seconds standard cine storage
- SSD-500G
- SSD-1T/2T(optional)
- HDD-1T/2T/4T (optional)
- USB Flash Drive(optional)
- Integrated DVDRW
- Integrated black/white thermal video printer slot
- 3D Mesh(optional)
- Patient information database
- Image archive on hard drive
- Quick store to USB memory stick
- Quick store to hard drive
- Quick print to B/W and color thermal

video printer

- Network storage and printing
- Full measurement and analysis

package

- Real time auto wave Doppler track and calculations
- Vascular calculations
- Cardiac calculations
- OB calculations and tables
- Gynecological calculations
- Urological calculations
- Renal calculations
- Volume calculations
- 3D Volume calculations(optional)
- Wireless networking for easy data sharing, storage and printing
- Bluetooth for image data transfer
- Gel Warmer(optional)
- Image data transfer directly by E-Mail with network access
- Up-to-date connectivity and data management solutions, wireless, LAN, Bluetooth, E-Mail, integrated database
- DICOM compatibility
- Four active probe ports
- 7 USB ports
- 8 TGC slides
- Average 4 multiple adjustable frequency in every probe and mode
- Up to 512 line density
- 1 DVI-D interface
- 1 Audio in interface, 1 Audio out interface
- 1 Speaker interface
- 1 RJ45 interface

1.5 Language support

- Software: Chinese, English, German, Greek, Malay, Portuguese, Romanian, Spanish, Swedish, French, Polish, Russian, Uyghur, Ukraine, Italian, Czech, Hungarian
- Keyboard input: Chinese, English, German, Greek, Malay, Portuguese, Romanian, Spanish, Swedish, Polish, Norwegian, Danish, Finnish, French, Russian, Italian, Czech, Cambodia, Polski, Ukraine
- Control panel overlay: English
- User manual: Chinese, English, German, Russian, Portuguese, Spanish, Italian, French

2. Ergonomics

- Unique human oriented design for comfort and convenience
- Fully articulating 23.8 inch high resolution flat panel display
- Easy access DVD media drive
- 4 easy access transducer ports
- 5 transducer holders (removable for easy cleaning, include one endocavity holder)
- Integrated touchable alphabetic keyboard
- Simple, easy and effective cable

management structure

2.1 Keyboard

- Intuitive, configurable and touchable interactive operation interface
- Ergonomic hard keys for general ultrasound operations
- 8 TGC slides, functionality at any depth
- Backlight keys
- Keyboard adjustable
 - Swivel range: -45° - 45°
 - Down/up range: 150 mm

2.2 Touch screen

Highly sensitive 15.6 inch LED technology touch screen

- Resolution: 1920*1080 pixels
- Intuitive, configurable and touchable interactive operation interface
- Ergonomic hard keys for general ultrasound operations
- Touch Screen adjustable
 - Tilt range: 45° - 90°

2.3 Image display screen

- 23.8 inch high resolution LED technology, pixel resolution:

1920x1080

- Brightness, contrast and color temperature adjustment
- View angle : -178° ~ 178°
- Number of color: 16.7M
- Adjustable Gamma curve optimization for dedicated applications
- Multifunctional support arm design
- Independent tilt and swivel adjustment
 - Swivel range: ± 180 degrees
 - Tilt range: $-20-90$ degrees
 - Up/down : 80mm

2.4 Wheels

- Diameter: 125mm
- Front castor (2 ea): Total lock
- Rear castor (2 ea): Total lock

2.5 Touch gestures

- Swipe down/up: display/remove projected image on touch screen
- Swipe horizontally: page up/down or review images/cine loops one by one
- Swipe from left edge to right: display hidden menu on projected image.
- Image parameter adjustment
- Measurement on projected image on touch screen
- Rotate or erase on projected 3D/4D image on touch screen

2.6 System boot-up

- Boot-up from shut-down: about 80sec
- Shut-down: about 9 sec

2.7 Comments

- Supports text input and arrow
- Support freehand marking on touch screen
- Adjustable text size and arrow size
- Supports home position
- Covers various application
- User customizable

2.8 Bodymark

- More than 215 bodymarks for versatile application
- User customizable

2.9 Peripherals

- B&W thermal video printer: Sony UP-D898MD(optional)
- Color thermal video printer: Sony UP-D25MD (optional)

2.10 Dimensions and Weight

- Height: 1420mm
- Width: 605mm
- Depth: 940mm
- Net Weight: 60kg

2.11 Electrical Power

- Voltage: 100-240V AC
- Frequency: 50/60Hz
- Power: < 700VA for console only
- Support built in battery(optional)
 - Scan time in B Mode: about 1h
 - Charging time: ≤4h

2.12 Operating Environment

- Ambient temperature: 10-40° C
- Relative humidity: 30-75%
- Atmospheric pressure: 700hPa-1060hPa

2.13 Storage & Transportation

Environment

- Ambient temperature: -10-50° C
- Relative humidity: 10%-90% (no condensation)
- Atmospheric pressure: 700hPa-1060hPa

3. Transducers

3.1 Transducer Technology

- Xcen technology for wideband frequency
- Pure wave technology for high resolution imaging
- Unique and high technical Xcen probe connector to adapt all different type of VINNO product models

3.2 Transducer types

- Convex array
- Linear array
- Phase array
- 4D probe
- Endocavity probe

3.3 Transducer selection

- Electronic switching of transducers
- User customizable imaging presets for each transducer and application
- Automatic dynamic receiving focus in all transducers
- Multiple adjustable transmit focal zone, up to 4 focal zone

F2-5C Broadband Curved Array

- Application: Abdomen, OB/Gyn, Urology, Pediatric
- Transducer Elements:128
- Physical Footprint: 72mm × 27mm
- Footprint: 17mm × 64mm
- Convex Radius: 60mm
- Field Of View: 59degree
- B Mode Frequency : 2.0-5.5MHz
- Harmonic Mode Frequency: 1.0-8.0MHz
- Center Frequency: 3.4MHz
- CF Mode Frequency: 2.0-4.0MHz
- PW Mode Frequency: 2.0-4.0MHz
- Pulsed Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in 2D, Harmonic, Color Doppler and Wave Doppler Modes
- Reusable Biopsy Guide available

S2-9C Broadband Curved Array

- Application: Abdomen, OB/Gyn, Urology, Pediatric
- Transducer Elements:128
- Physical Footprint: 88mm × 31mm
- Footprint: 18mm × 67mm
- Convex Radius: 60mm
- Field Of View: 66degree
- B Mode Frequency: 2-5.5MHz
- Harmonic Mode Frequency: 1.0-8.0MHz
- Center Frequency: 3.5MHz
- CF Mode Frequency: 1.8-4.0MHz
- PW Mode Frequency: 1.8-4.0MHz

- Pulsed Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in 2D, Harmonic, Color Doppler and Wave Doppler Modes
- Reusable Biopsy Guide available

X2-6C Broadband Curved Array

- Application: Abdomen, OB/Gyn, Urology, Pediatric
- Transducer Elements:192
- Physical Footprint: 76mm × 27mm
- Footprint: 16.8mm × 70mm
- Convex Radius: 60mm
- Field Of View: 65degree
- B Mode Frequency : 2.0-5.5MHz
- Harmonic Mode Frequency: 1.5-6.0MHz
- Center Frequency: 4.0MHz
- CF Mode Frequency: 1.8-4.0MHz
- PW Mode Frequency: 1.8-4.0MHz
- Pulsed Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in 2D, Harmonic, Color Doppler and Wave Doppler Modes
- Reusable Biopsy Guide available

D2-6C broadband curved array volume probe

- Application: Abdomen, OB/Gyn, Urology
- Transducer Elements:128
- Physical Footprint: 75.5mm × 49.2mm

- Convex Radius: 40mm
- Field Of View: 75degree
- B Mode Frequency : 3.0-5.5MHz
- Harmonic Mode Frequency: 3.0-6.0MHz
- Center Frequency: 4.0MHz
- CF Mode Frequency: 2.5-4.0MHz
- PW Mode Frequency: 2.5-4.0MHz
- Pulsed Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode, 3D/4D Grayscale and 3D Color Modes
- Multi-Imaging Frequency Setting in 2D, 3D/4D, Harmonic, Color Doppler and Wave Doppler Modes

G3-9M broadband micro convex array

- Application: Pediatric, Abdomen, Cardiac
- Transducer Elements:128
- Physical Footprint: 34.2mm × 28.7mm
- Footprint: 11.2mm × 25mm
- Convex Radius: 15mm
- Field Of View: 103degree
- B Mode Frequency : 5.0-10.0MHz
- Harmonic Mode Frequency: 4.0-12.0MHz
- Center Frequency: 6.5MHz
- CF Mode Frequency: 4.0-5.0MHz
- PW Mode Frequency: 4.0-5.0MHz
- Pulsed Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in 2D, Harmonic, Color Doppler and Wave Doppler Modes

- Reusable Biopsy Guide available

F4-9E broadband micro convex endocavity array

- Application: OB/Gyn, Urology
- Transducer Elements:128
- Physical Footprint: 32.5mm x 44.2mm
- Footprint: 10.7mm × 21mm
- Convex Radius: 10mm
- Field Of View: 150degree
- B Mode Frequency: 5.0-10.0MHz
- Harmonic Mode Frequency: 4.0-11.0MHz
- Center Frequency: 6.5MHz
- CF Mode Frequency: 4.0-5.0MHz
- PW Mode Frequency: 4.0-5.0MHz
- Pulsed Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in 2D, Harmonic, Color Doppler and Wave Doppler Modes
- Reusable Biopsy Guide available

G4-9E broadband micro convex endocavity array

- Application: OB/Gyn, Urology
- Transducer Elements: 128
- Physical Footprint: 40mm x 32.4mm
- Footprint: 9mm × 18.1mm
- Convex Radius: 11.5mm
- Field Of View: 136degree
- B Mode Frequency: 5.0-10.0MHz
- Harmonic Mode Frequency: 4.0-11.0MHz

- Center Frequency: 6.5MHz
- CF Mode Frequency: 4.0-5.0MHz
- PW Mode Frequency: 4.0-5.0MHz
- Pulsed Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in 2D, Harmonic, Color Doppler and Wave Doppler Modes
- Reusable Biopsy Guide available

G4-9EV broadband micro convex endocavity array

- Application: OB/Gyn, Urology
- Transducer Elements: 160
- Physical Footprint: 65mm x 28mm
- Footprint: 11mm × 25mm
- Convex Radius: 12mm
- Field Of View: 136degree
- B Mode Frequency: 5.0-10.0MHz
- Harmonic Mode Frequency: 4.0-11.0MHz
- Center Frequency: 6.5MHz
- CF Mode Frequency: 4.0-5.0MHz
- PW Mode Frequency: 4.0-5.0MHz
- Pulsed Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in 2D, Harmonic, Color Doppler and Wave Doppler Modes
- Reusable Biopsy Guide available

X4-9E broadband micro convex endocavity array(crank and straight handle)

- Application: OB/Gyn, Urology

- Transducer Elements:192
- Physical Footprint: 32.4mm x 40mm
- Footprint: 10mm × 18mm
- Convex Radius: 8.8mm
- Field Of View: 180degree
- B Mode Frequency: 5.0-10.0MHz
- Harmonic Mode Frequency:
4.0-11.0MHz
- Center Frequency: 6.5MHz
- CF Mode Frequency: 4.0-5.0MHz
- PW Mode Frequency: 4.0-5.0MHz
- Pulsed Wave Doppler, Color Doppler,
Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in
2D, Harmonic, Color Doppler and Wave
Doppler Modes
- Reusable Biopsy Guide available

D4-9E broadband micro convex 4D endocavity array

- Application: Ob/Gyn, Urology
- Transducer Elements:148
- Physical Footprint:32.5mm x
39.5mm
- Convex Radius: 10mm
- Field Of View: 141degree
- B Mode Frequency: 5.0 - 10.0MHz
- Harmonic Mode Frequency:
4.0-11.0MHz
- Center Frequency: 6.5MHz
- CF Mode Frequency: 4.0-5.0MHz
- PW Mode Frequency: 4.0-5.0MHz
- Pulsed Wave Doppler, Color Doppler,
Power Doppler, Harmonic, B-Mode,
3D/4D Grayscale

- Multi-Imaging Frequency Setting in
2D, 3D/4D, Harmonic, Color Doppler
and Wave Doppler Modes

BP4-9 Biplane array

BP4-9L

- Application: Gyn, Urology
 - Transducer Elements:128
 - Physical Footprint: 30mm x 26mm
 - Footprint: 8mm × 60mm
 - Aperture Size: 26mm
 - B Mode Frequency Range:
6.0-12.0MHz
 - Harmonic Mode Frequency:
6.0-15.0MHz
 - Center Frequency: 7.3MHz
 - CF Mode Frequency: 3.0-6.3MHz
 - PW Mode Frequency: 3.0-6.3MHz
 - Pulsed Wave Doppler, Color Doppler,
Power Doppler, Harmonic, B-Mode
- ##### BP4-9C

- Application: Gyn, Urology
- Transducer Elements:128
- Physical Footprint: 30mm x 26mm
- Footprint: 11mm × 24mm
- Convex Radius: 10mm
- Field Of View: 150degree
- B Mode Frequency: 5.0-10.0MHz
- Harmonic Mode Frequency:
4.0-9.0MHz
- Center Frequency: 6.5MHz
- CF Mode Frequency: 4.0-5.0MHz
- PW Mode Frequency: 4.0-5.0MHz
- Pulsed Wave Doppler, Color Doppler,
Power Doppler, Harmonic, B-Mode

- Reusable Biopsy Guide available

F4-12L broadband linear array

- Applications: Vascular, Small Parts
- Transducer Elements:128
- Physical Footprint: 52.5mm × 25mm
- Footprint: 9mm × 44mm
- Aperture Size: 38.4mm
- B Mode Frequency : 6.0 -12.0MHz
- Harmonic Mode Frequency:
6.0-18.0MHz
- Center Frequency: 7.3MHz
- CF Mode Frequency: 3.0-6.3MHz
- PW Mode Frequency: 3.0-6.3MHz
- Pulsed Wave Doppler, Color Doppler,
Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in
2D, Harmonic, Color Doppler and Wave
Doppler Modes
- Reusable Biopsy Guide available

X4-12L broadband linear array

- Applications: Vascular, Small Parts
- Transducer Elements:192
- Physical Footprint: 53.2mm ×
23.4mm
- Footprint: 6.7mm × 40mm
- Aperture Size: 38.4mm
- B Mode Frequency : 6.0 -12.0MHz
- Harmonic Mode Frequency:
6.0-17.0MHz
- Center Frequency: 7.3MHz
- CF Mode Frequency: 3.0-6.3MHz
- PW Mode Frequency: 3.0-6.3MHz
- Pulsed Wave Doppler, Color Doppler,

Power Doppler, Harmonic, B-Mode

- Multi-Imaging Frequency Setting in
2D, Harmonic, Color Doppler and Wave
Doppler Modes

- Reusable Biopsy Guide available

X6-16L broadband linear array

- Applications: Vascular, Small Parts
- Transducer Elements:192
- Physical Footprint: 52.8mm × 26mm
- Footprint: 6.7mm × 40mm
- Aperture Size: 38.4mm
- B Mode Frequency: 7.3 -14.0MHz
- Harmonic Mode Frequency:
8.0-19.0MHz
- Center Frequency: 10.0MHz
- CF Mode Frequency: 5.0-13.0MHz
- PW Mode Frequency: 5.0-13.0MHz
- Pulsed Wave Doppler, Color Doppler,
Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in
2D, Harmonic, Color Doppler and Wave
Doppler Modes
- Reusable Biopsy Guide available

X3-10L Low-frequency linear array

- Applications: Musculoskeletal,
Peripheral Vascular
- Transducer Elements:192
- Physical Footprint: 61mm × 26mm
- Footprint: 11mm × 50mm
- Aperture Size: 46mm
- B Mode Frequency: 4.0-10.0MHz
- Harmonic Mode Frequency:
3.0-15.0MHz

- Center Frequency: 6.0MHz
- CF Mode Frequency: 3.0-6.3MHz
- PW Mode Frequency: 3.0-6.3MHz
- Pulsed Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in 2D, Harmonic, Color Doppler and Wave Doppler Modes
- Reusable Biopsy Guide available

U5-15LE broadband linear array

- Applications: Small Parts, Breast, Vascular
- Transducer Elements:256
- Physical Footprint: 99mm×26mm
- Footprint: 7.5mm×55mm
- Aperture Size: 51mm
- B Mode Frequency : 6.0 -12.0MHz
- Harmonic Mode Frequency: 8.0-15.0MHz
- Center Frequency: 8.5MHz
- CF Mode Frequency: 5.0-8.0MHz
- PW Mode Frequency: 5.0-8.3MHz
- Pulsed Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in 2D, Harmonic, Color Doppler and Wave Doppler Modes
- Reusable Biopsy Guide available

I7-18L broadband linear array (Hock Stick)

- Applications: Vascular, Small Parts
- Transducer Elements:128
- Physical Footprint: 31.8mm x 10mm

- Footprint: 4.9mm×28mm
- Aperture Size: 25.6mm
- B Mode Frequency: 6.0-12.0MHz
- Harmonic Mode Frequency: 8.0-20.0MHz
- Center Frequency: 8.5MHz
- CF Mode Frequency: 5.0-13.0MHz
- PW Mode Frequency: 5.0-13.0MHz
- Pulsed Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in 2D, Harmonic, Color Doppler and Wave Doppler Modes

X9-22L broadband linear array

- Applications: MSK, Nerve, Small Parts
- Transducer Elements:192
- Physical Footprint: 45mm x 25.4mm
- Footprint: 5.1mm×32mm
- Aperture size: 28.8mm
- B mode Frequency : 9.0-18.0MHz
- Harmonic Mode Frequency: 12.0-23.0MHz
- Center Frequency: 15.0MHz
- CF Mode Frequency: 8.3-16.7MHz
- PW Mode Frequency: 8.3-16.7MHz
- Pulsed Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in 2D, Harmonic, Color Doppler and Wave Doppler Modes

X10-23L broadband linear array

- Applications: Small Parts

- Transducer Elements:128
- Physical Footprint: 49.6mm x 22.2mm
- Footprint: 4mm × 16mm
- Aperture Size: 12.8mm
- B Mode Frequency : 10.0-20.0MHz
- Harmonic Mode Frequency: 15.0-24.0MHz
- Center Frequency: 13.0MHz
- CF Mode Frequency: 11.1-16.7MHz
- PW Mode Frequency: 11.1-16.7MHz
- Pulsed Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in 2D, Harmonic, Color Doppler and Wave Doppler Modes

G1-4P phased array

- Applications: Cardiac, Abdomen, TCD
- Transducer Elements: 64
- Physical Footprint: 34.2mm × 28.7mm
- Footprint: 15mm × 22mm
- Aperture Size: 18mm
- Field Of View: 90degree
- B Mode Frequency : 2.0-3.5 MHz
- Harmonic Mode Frequency: 1.0-6.0MHz
- Center Frequency: 2.8MHz
- CF Mode Frequency: 1.7-3.3MHz
- PW Mode Frequency: 1.7-3.3MHz
- Pulsed Wave Doppler, Continuous Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in

2D, Harmonic, Color Doppler and Wave Doppler Modes

- Reusable Biopsy Guide available

S1-6P phased array

- Single Crystal Technology
- Applications: Cardiac, Abdomen, TCD
- Transducer Elements: 96
- Physical Footprint: 36mm x 29mm
- Footprint: 16mm × 23mm
- Aperture Size: 19.2mm
- Field Of View: 90degree
- B Mode Frequency : 2.0-5.0 MHz
- Harmonic Mode Frequency: 1.0-8.0MHz
- Center Frequency: 3.5MHz
- CF Mode Frequency: 1.7-3.3MHz
- PW Mode Frequency: 1.7-3.3MHz
- Pulsed Wave Doppler, Continuous Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in 2D, Harmonic, Color Doppler and Wave Doppler Modes
- Reusable Biopsy Guide available

G3-10P phased array

- Application: Pediatric Cardiology, Abdomen
- Transducer Elements: 64
- Physical Footprint: 33mm x 33mm
- Footprint: 12mm × 18.6mm
- Aperture Size: 10.2mm
- Field Of View: 90degree
- B Mode Frequency: 3.0-8.0 MHz

- Harmonic Mode Frequency: 3.0-10.0MHz
- Center Frequency: 5.0MHz
- CF Mode Frequency: 3.3-5.7MHz
- PW Mode Frequency: 3.3-5.7MHz
- Pulsed Wave Doppler, Continuous Wave Doppler, Color Doppler, Power Doppler, Harmonic, B-Mode
- Multi-Imaging Frequency Setting in 2D, Harmonic, Color Doppler and Wave Doppler Modes

G1-3R Pencil Probe

- Applications: TCD
- Transducer Elements:2
- Physical Footprint: 17.1mm × 17.1mm

4. Advanced Imaging controls

4.1 VFusion

- Available on all transducers and for 2D, 3D/4D (except phase array probe)
- Operate in conjunction with VSpeckle, harmonic imaging

4.2 VSpeckle

- Available on all transducers and for 2D, 3D/4D

- Virtually eliminate speckle noise artifact and dynamically enhances tissue margins
- Selectable multiple levels of speckle noise reduction and smoothing
- Operates in conjunction with VFusion and harmonic imaging

4.3 VTissue

- Advanced imaging processing to adapt to the speed of the ultrasound variation in different tissue
- Improved detail resolution and conspicuity of lesions
- Presentable sound and speed in different applications
- One touch operation to ease diagnosis
- Better detection in diffuse lesions of organs

4.4 3D/4D

4.4.1 3D/4D HQ Grad

- Amazing high image quality
- Extreme realistic rendering images
- Similar operation as normal Rendering

4.4.2 Inversion mode

- This render mode is used to

display anechoic structures such as vessels

- It invert the gray values of the rendered image, such as black image information become white and vice versa

4.4.3 Magic Cut(optional)

- Ability to edit images, make possible to cut away structure obstructing the view in the ROI
- Several cutting methods available
- Have quality index to indicate if there is proper external force

4.4.4 Niche view

- Display 3 orthogonal planes centered on ROI
- Use Depth to translate the selected plane
- Each imaging plane or Niche image can be selected using image reference

4.4.5 Free View(optional)

- Provide any plane view to visualize the internal tissue information
- Improve the contrast resolution to facilitate the detection of diffuse lesions in organs

4.4.6 3D Smart Face(optional)

An intelligent tool for fetal face optimization. This tool detects the fluid/tissue interface and smartly removes noise in front of the baby inside the ROI, to obtain an optimal baby face.

- Use Auto key on the keyboard to remove the obstacle in front of the baby
- Only works on 3D Render
- Can not use this feature together with MagicCut

4.4.7 VNavIn (optional)

A tool that navigates inside the 3D volume data and projects an inside-out perspective image that displays the inner most structures like virtual endoscopy

- This feature is useful in body structures which are surrounded by fluid, like gynecology, obstetrics, abdomen, vascular, or any other fluid-filled areas.
- Two ways to perform VNavIn Auto and Manual
- Depth mapping to enhance the depth perception

4.5 Tissue Doppler (TD)

- Present wall motion spectrum by

using Doppler principle

- Provide wall motion direction and velocity information

4.6 Tissue Velocity Imaging (TVI)

- Color codes the velocities in tissue
- Present tissue color imaging by using Doppler principle
- This color image is overlaid onto the 2D image
- Captures low flow but high amplitude signals associated with wall motion

4.7 Stress Echo(optional)

- Stress echo is a non-invasive, dynamic evaluation of myocardial structure and its function under an external stress(exercise or pharmacology)
- 12 Ready to use templates (max 8 stages * 6 views) Editable
- User definable template
- Re-arrange & Select default template
- 10 View names available
- 14 Stage names are available (can add user defined stage name)
- One Touch Shuffle (Stage / View)
- Touch & Compare any view of stage
- Systole only review

4.8 Strain Imaging(optional)

- Auto-ROI (after selecting Mitral Valve Plane)
- Adjust Segment-wise (Longitudinal strain)
- Adjust Segment-wise and Rotate whole ROI (Radial & Circumf. Strain)
- ECG to select heart cycle
- View based Bulls Eye view
- Result type (Peak Strain or Peak Time)Parameter type (L Strain & C Strain)

4.9 Tissue Velocity M mode(TVM)

(optional)

- Color codes the velocities in tissue
- Present wall motion spectrum based on tissue moving
- This color image is overlaid onto the 2D image
- Captures low flow but high amplitude signals associated with wall motion

4.10 Multi-angle M mode

- Sample on moving tissue from multi-angle
- Present wall motion spectrum based on tissue moving

4.11 Curved M mode

Draw the route of the sample line freely and obtain the corresponding anatomical M-mode. This might be helpful to obtain myocardial wall motion.

- Color curved M mode is available
- TVI & M mode is available

4.12 VAid (optional)

- An artificial intelligent detection tool for breast/Thyroid imaging
- Works in real-time detection, as well as on the stored (single or cine) imagele

5. Imaging modes

5.1 2D Imaging

- Pre-defined ATGC (adaptive temporal gain compensation) curves optimized for consistently excellent imaging
- Display format: Single, Dual, Quad
- B/M acoustic output: 10-100%
- Select between 0 to 4 transmit

focal zones

- Reverse function: on/off
- 2D optimization: on/off
- Centerline: on/off
- L/R flip and U/D flip: on/off
- VFusion : \cong 7steps
- VSpeckle : \cong 6 steps
- Harmonic imaging both tissue harmonic and phase inversion
- Cineloop image review
- Selectable 2D line density
- Dual imaging with independent cineloop
- 256(8 bit) gray level
- FULL screen imaging to larger image size
- Multi frequency: \cong 4 levels, probe dependent
- Gray filter: \cong 7 steps
- Persistence: \cong 8steps
- Selectable image angles, probe Dependent
- Gain: 0-100%
- Dynamic range: 30-280 db
- VSharpen(enhance edge contrast) : \cong 8steps
- Smooth(improve spatial resolution): \cong 11steps
- EdgeEnhance (improve detail information and contrast): \cong 6steps
- Gray Map: \cong 32types
- Tint Map: \cong 24types
- TGC: 8 slides on control pannel
- SGC: 8 ponds on touch pannel
- TI heat index: TIB, TIS, TIC

- Rotation:
0° ,90° ,180° ,270°
- Zoom(up to 10×)

5.2 Harmonic Imaging

- Supports both tissue harmonic and phase inversion imaging (transducer and frequency dependence)
- Second harmonic processing to reduce artifacts and improve image clarity
- Maximize detail resolution and enhance contrast
- Available on all imaging transducers
- Extends high performance imaging capabilities to all patient body types

5.3 M mode

- Selectable sweeping rates, 10steps
- Selectable display format prospective or retrospective (V2/3, V1/3, V1/2, H1/2, H3/4, full screen)
- Chroma colorization with multiple color maps
- Cineloop review for retrospective analysis of M-mode data
- 256 gray levels
- Acoustic output: 10%-100%
- Gray filter: \cong 7steps

- Dynamic range: 108db-128db, 2db/step
- Vsharpen: \cong 6steps
- Gray Map: \cong 32types
- Tint Map: \cong 24types
- Gain: 0-100%
- Color M mode: available
- MultiAngle: available

5.4 Color Doppler mode

- Available on all imaging transducers
- Automatically adapts transmit and receive bandwidth processing based on the color box position
- Cineloop review with full playback control
- Color flow M mode display for tissue motion and flow velocity
- Revert: on/off
- Selectable baseline, line density, flash reduction, persistence, maps, frequency, PRF, wall filter, packet size, color level, sensitivity, focus position, acoustic power, and smooth
- FULL screen imaging to larger image size
- L/R flip and U/D flip: on/off
- Frequency: \cong 4steps, depend on probes
- Baseline: 0-100%
- Acoustic power: 5% -100%
- Line density: \cong 5 steps

- Flash reduction: \cong 5 steps
- Persistence: \cong 7 steps
- Color Map: \cong 33types
- Smooth : \cong 7steps
- Sensitivity: \cong 5 steps
- Transparency: \cong 6steps
- Color level: \cong 15steps
- Packet size: \cong 7 steps
- Color gain: 0-100%
- Adjustable region of interest
- Region of interest
- Baseline invert
- Simultaneous mode during PW mode
- Zoom

5.5 Power Doppler mode

- High sensitive mode for small vessel visualization
- Available on all transducers
- Cineloop review
- Display format: Single, Dual, Quad
- Selectable line density, flash reduction, persistence, maps, frequency, PRF, wall filter, packet size, color level, sensitivity, focus position, acoustic power, and smooth
- Color maps: \cong 24 types
- Color levels: \cong 11 steps
- Sensitivity: \cong 5steps
- Smooth: \cong 7steps
- Persistence: \cong 7steps

- Individual controls for gain
- Adjustable region of interest

5.6 Pulsed Wave (PW) Doppler

- Angle correction with automatic velocity scale adjustment
- Normal, invert display around horizontal zero line
- Auto optimization: on/off
- Invert: on/off
- Selectable display format prospective or retrospective (V2/3, V1/3, V1/2, H1/2, H3/4, full screen)
- Selectable gray filter, dynamic range, frequency, PRF, wall filter, baseline, angel correct, sample volume
- Gray filter: \cong 6steps
- Dynamic range: 108db-128db
- Baseline: 5%-95%
- Sample volume: 0mm-28mm
- Angle correct: $-80^{\circ} \sim 80^{\circ}$
- Trace Sensitivity: \cong 21steps
- Audio Volume: \cong 27steps
- Spectrum Optimize: \cong 28steps
- Gray map: \cong 13types
- Tint map: \cong 11types
- Selectable sweep speeds: \cong 10 steps
- Maximum velocity range: 12m/s
- PW acoustic output: 5%-100%
- Trace direction: above, below, above and below

- Trace type: max, mean, max and mean
- Cardiac cycle: 1-5
- Selectable low frequency signal filtering with adjustable wall filter settings
- Selectable chroma colorization maps
- Auto function to optimize spectral Doppler display
- Digitally enhanced stereo output
- 256 gray levels
- Post-processing in frozen mode includes map, baseline, invert and chroma
- Simultaneous or duplex mode of operation
- Simultaneous 2D, color Doppler, pulsed Doppler
- High PRF capability in all modes including duplex and triplex

5.7 Continuous Wave Doppler (CW)

- User can measure distance and area
- Cineloop review
- Selectable gray filter, dynamic range, PRF, wall filter, baseline, angle correct, sample volume, acoustic power, etc.
- Gray filter: \cong 6 steps
- Dynamic range: 108db-128db
- Angle correct: $-80^{\circ}\sim 80^{\circ}$
- Trace Sensitivity: \cong 21 steps
- Audio Volume: \cong 27 steps

- Spectrum Optimize: \cong 28 steps
- Gray map: \cong 13 types
- Tint map: \cong 11 types
- acoustic output: 5%-100%

5.8 Elastography imaging

- Shows the spatial distribution of tissue elasticity properties in a region of interest to estimate the strain before and after tissue distortion caused by external force
 - The strain estimation is scaled by color to have smooth distribution display
 - Have quality index to indicate if there is proper external force
- Precision: 0,1,2,3,4
- Resolution: 0,1,2,3,4
- Sensitivity: 0-10
- Transparency: \cong 13steps
- Smooth: \cong 7steps
- Line density: \cong 7steps
- Persistence: \cong 7steps
- Map: EIO
- Display format: Single, Dual, Quad

5.9 3D/4D

- 3D/4D rotation
- Grayscale imaging controls
- Selectable rendering

Approaches: HQ Surface, HQ Grad, HQ Silhouette, Surf Texture, Surf Smooth, Grad Light, Surf HDR, Trans Max, X-ray, Transp Min, Light

- Unique high quality rendering algorithm
- Review volume
- Volume Angle: 15%-85%
- Quality: low, mid, good, high, best
- Threshold: 0-255 256Steps
- Transparency: 0.1-2, 0.1/step
- Display format: single, dual, MRP, Quad
- Image Reference: A, B, C, 3D
- Flip: 0° , 90° , 180° , 270°
- View: Front/Back, Back/Front; Left/Right, Right/Left; Up/Down, Down/Up
- Rotation Direction: X, Y, Z
- 3D Map: \cong 10 types
- Tint maps: \cong 24Types
- Gray maps: \cong 32Types
- 2D VSpeckle: \cong 4types
- 3D VSpeckle: \cong 4types
- Render Type: Gray, GrayInv

MCUT

- Slice Number: 2x2, 3x3, 4x4, 5x5
- Max Slice Number: 25
- Gray Map: \cong 32 types
- Tint Map: \cong 24 types
- Cut plane: A, B, C
- Rotation Direction: X, Y, Z
- Volume Angle: 15-85°
- Interval: 1mm-20mm,

0.5mm/step

- Quality: low, mid, good, high, best

Free view(optional)

- Direction: X, Y, Z
- Route: curve, straight line
- Reference image: A, B, C
- Slice thickness: 0mm-20mm
- Active line: 1, 2, 3
- Mix: 10-90
- Threshold: 256 steps
- Transparency: 0.1-2.0, 0.1/step

Magic cut(optional)

- Erase mode: Inside Lasso, Outside Lasso, big circle, small circle
- Erase type: trace, rectangle, ellipse
- Rotation direction: X, Y, Z

Niche view

- Model type: upper, lower
- Display format: single, quad
- Rotation direction: X, Y, Z
- Image reference: A, B, C, N

5.10 PView (optional)

- Real time extended field of view composite imaging
- Ability to back up and realign the image during acquisition
- Full zoom, cinelooop review and image rotation capabilities

5.11 TView

- Expand view of scanning

5.12 Auto optimization

- Intelligent one button automatic optimization in 2D and Doppler modes
- Automatically adjust PRF and baseline in Doppler

5.13 Live Track (including angle free flow)(optional)

- Once the function is enabled, system automatically tracks the vascular and adjust the ROI position、angle in CF and PW mode, the sample volume、steer、correct angle
- Applications: carotid
- Probe:
X4-12L,X6-16L,X3-10L,U5-15LE,F4-12L

5.14 Tissue Doppler Imaging (TD)

- Present wall motion spectrum by using Doppler principle
- Provide wall motion direction and velocity information
- Available on all sector transducer for cardiac imaging
- Selectable frequency, PRF, wall filter

- Gain
- Sweep speed: $\cong 10$ steps
- Baseline: 5%-95%
- Angle correct: $\pm 80^\circ$
- Sample volume: 0.5mm-10mm
- Spectrum optimize: $\cong 20$ steps
- Acoustic power: 5%-100%
- Dynamic range: 92db-128db
- Trace sensitive: $\cong 21$ steps
- Gray filter: $\cong 6$ steps
- Audio volume: $\cong 27$ steps
- Mode: max, mean, max and mean
- Direction: above, below, above and below
- Heart cycle: 1-5
- Gray map: $\cong 13$ types
- Tint map: $\cong 11$ types

5.15 Tissue Velocity Imaging (TVI)

- Color codes the velocities in tissue
- Present tissue color imaging by using Doppler principle
- This color image is overlaid onto the 2D image
- Captures low flow but high amplitude signals associated with wall motion
- Available on all sector transducer for cardiac imaging
- Tissue velocity M mode display for wall motion(optional)

- Gain
- Velocity
- Color level: \cong 16steps
- Transparency: \cong 13steps
- Smooth: \cong 7steps
- Line density: \cong 3steps
- Persistence: \cong 7steps
- Color map: \cong 10types

5.16 Tissue Velocity M mode (TVM)

(optional)

- Color codes the velocities in tissue
- Present wall motion spectrum based on tissue moving
- This color image is overlaid onto the 2D image
- Captures low flow but high amplitude signals associated with wall motion
- Selectable frequency, PRF, Focus position
- Baseline: 5%-95%(M), 0%-100%(CF)
- Color level: \cong 16steps
- Transparency: \cong 13steps
- Packet size: 3,4,5,6
- Acoustic power: 5%-100%(CF), 10%-100%(M)
- Display format: Single, Dual, Quad

6. Touch Panel Interface

6.1 2D mode

- New patient
- BodyPattern
- Archive
- Comments
- End exam
- Report
- Sys setting
- Probe&App
- Pview
- Tview
- Fullscreen
- L/R
- U/D
- Center line
- VTissue
- VSpeckle
- VFusion
- Gray Filter
- Persistence
- Display Format
- Image reference
- Maps
- Frequency
- Focus position
- Focus#
- Dynamic Range
- Line density
- VSharpen
- Biopsy
- Image angle

- Focus width
- Smooth
- Acoustic power
- Elastosonography
- EdgeEnhance
- NeedleEnhance
- SGC

6.2 M Mode

- New patient
- BodyPattern
- Archive
- Comments
- End exam
- Report
- Sys setting
- Probe&App
- Maps
- Dynamic range
- Acoustic power
- Sweep speed
- Gray filter
- VSharpen
- ECG

6.3 CF mode

- New patient
- BodyPattern
- Archive
- Comments
- End exam
- Report
- Sys setting

- Probe&App
- Invert
- Full Screen
- L/R
- U/D
- Baseline
- Flash Reduction
- Line density
- Persistence
- Display format
- Sync display
- Transparency
- Image reference
- Maps
- Frequency
- PRF
- Wall filter
- Packet size
- Colorlevel
- Sensitivity
- Focus position
- Acoustic power
- Smooth

6.4 PW/CW mode

- New patient
- BodyPattern
- Archive
- Comments
- End exam
- Report
- Sys setting
- Probe&App

- Invert
- Duplex/Triplex
- Sweep speed
- Gray filter
- Dynamic range
- Trace sensitive
- Auto trace
- Mode/direction
- Maps
- Frequency
- PRF
- Wall filter
- Baseline
- Steer
- Sample volume
- Audio Volume
- Spectrum optimize
- Acoustic power

6.5 3D mode

- Comments
- Probe&App
- BodyPattern
- Probe&App
- Back to 2D
- Start3D
- Render
- Display format
- Image reference
- View
- Gray map
- VSpeckle
- Quality

- Threshold
- Transparency
- Volume angle
- Auto rotate (after data acquisition)
- Movement step (after data acquisition)
- HQ Grad(after data acquisition)
- Rotation angle (after data acquisition)
- Rotation direction (after data acquisition)
- 3DMcut(after data acquisition)
- Magic Cut (after data acquisition)
- Free View(after data acquisition)

6.6 4D mode

- Comments
- Probe&App
- Body Pattern
- Probe&App
- Back to 2D
- Start 4D
- Rotation direction
- Render
- Display format
- Image reference
- View
- Gray map
- VSpeckle
- Quality
- Threshold
- Transparency
- Volume angle

7. System Feature

7.1 Display modes

- Simultaneous capability
 - 2D/PW/CW
 - 2D/CF or PDI
 - 2D/M
 - Dual, 2D/2D
 - Dual, 2D/2D+CF or PDI
 - Dual, duplex and triplex
 - Duplex and Triplex mode
 - Quad display in 3D/4D application
 - 25 slice images display in 3D/4D application
- Independent dual 2D/PW or CW
- Timed based sweep update mode

7.2 Display annotation

- Institution/hospital name
- Date: 3 types selectable, Year-Month-Day, Day-Month-Year, Month-Day-Year
- Time: 2 types selectable, 24hours and 12 hours
- Operator identification

- Patient: First Name, Middle Name, Family Name
- Patient identification: 20 characters
- Gestational age from LMP/BBT/DOC/IVF/GA/Avg.US
- VINNO image symbol: Ginkgo leaf
- Power output index
 - MI: mechanical index
 - TIS: thermal index soft tissue
 - TIC: thermal index cranial (Bone)
 - TIB: thermal index bone
- Probe orientation marker: coincide with a probe orientation marking on the probe
- Gray/color bar
- Measurement result window
- Probe type
- Application name
- Image depth
- Imaging parameters by mode
 - 2D/M mode: acoustic power output, gain, frequency, frame rate, dynamic range
 - Color mode: color acoustic power output, color gain, color flow frequency, PRF, wall filter
 - PW/CW mode: Doppler acoustic power output, Doppler gain, Doppler frequency, PRF, wall filter, sample depth, SV(PW)
- Scanline Gain Compensation(SGC) with 8 slides

adjustment

- Focus zone marker
- Body pattern
- PW and CW scale markers:

time/speed

- M scale markers: time/depth,

time

- System measurement display
- System message display
- Biopsy guide line
- Heart rate

7.3 Simple User Operation Interface

- Simple user interface and easy workflow, allows one step on probe & application switch, and intuitive user parameter control

7.4 Cineloop

- Acquisition, storage in memory and display of up to 30000 frames, 1500 seconds long of 2D, color and PW/CW images for review
- Available to decide StartFrame and EndFrame
- Frame by frame manual cine loop review or auto playback with variable speed:
400%, 200%, 100%, 60%, 50%, 40%, 20%
- Frame compare: displays one cine in dual format and allows frame by frame

compare side by side

- Acquisition, storage and replay of Doppler audio

7.5 Quick save feature

- The system provides quick save function through USB stick, internal/external HDD, DVD during or after exam
- Configurable saving file format, VRD (VINNO Raw Data), DICOM, AVI, MP4, PNG

7.6 ECG (optional)

- One 3-lead ECG input
- Gain, sweep rate and display position controls
- Automatic heart rate calculation and display
- Fault condition display

7.7 Archive

- Patient data input which include patient ID, First Name, Middle Name, Family Name, Data of Birth, sex, Perf.Physician, Ref.Physician, exam operator
- Physical data such as weight, height
- Patient exam management
- Patient exam images storage and

management

- Import VRD format data into the system from outside media, such as USB stick, external HDD, DVD
- Export patient data into outside medias

7.8 Report

- Automatically pull patient data into the report
- Automatically load measurement worksheet into the report
- Pull related exams' images into the report
- Write comments in the report
- Print report through network or local printer

7.9 Connectivity

- Standard connectivity features
 - Local print to on-board or off-board video printers through USB port
 - Page report print
 - Image export to removable media (DVD, external HDD, USB stick)
- Ethernet Network Connection
 - Cable connection
 - Wireless connection: need wireless routing adaptor
- Network linkage
 - Image export to

network storage servers

- DICOM export and retrieve (optional)
- Mobile data transfer solution by
 - Blue tooth
 - email
 - Hot point connection
- VCloud
- export 3D data for 3D printer(optional)
- Integrated DVDRW
 - Support standard DVD media
 - Data storage formats include VRD, DICOM, AVI, MP4, PNG
 - JPEG,BMP,PNG,VRD and DICOM images stored in disc can be recalled on the VINNO system
 - PNG and AVI/MP4 images can be played on normal computers
- On-board patient exam storage
- Direct digital storage of static image or cineloop images to internal hard disk drives
- Fully integrated user interface

7.10 Probes/application

- Selectable multiple applications
- Edit exist application preset
- Edit user defined preset
- Rename preset
- Return to factory preset
- Quick save user defined

parameters in related application

7.11 Safety Conformance

- Regulatory Notice:

This device is tested to meet all applicable requirements in relevant. According to Regulation (EU) 2017/745 concerning medical devices.

- Conformity to Standards:

- IEC

60601-1:2005/A1:2012+A2:2020

Medical electrical equipment - Part 1: General requirements for basic safety and essential performance

- IEC 60601-1-2:2014/A1:2020

Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests

- IEC

60601-1-6:2010/A1:2013+A2:2020

Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability

- IEC 60601-2-37:2007/A1:2015

Medical electrical equipment - Part 2-37: Particular requirements for the basic safety and essential performance of ultrasonic medical diagnostic and monitoring equipment

- IEC 61157:2007/A1:2013 Standard

means for the reporting of the acoustic output of medical diagnostic ultrasonic equipment

- ISO 10993-1:2018 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process

- IEC 62304:2006/A1:2015 Medical device software - Software life-cycle processes

- IEC 62366-1:2015/A1:2020 Medical devices - Application of usability engineering to medical devices

- WEEE according to 2012/19/EU

- RoHS according to 2011/65/EU

8. Measurement and Analysis

8.1 Measurement in different modes

8.1.1 Generic Measurement in 2D mode

- Depth
- Distance
- Perimeter
 - Length and width

method

- Ellipse method
- Polygon method
- Spline method
- Tracing method

- Area
 - Length and width method
 - Ellipse method
 - Polygon method
 - Spline method
 - Tracing method
- Volume
 - Single line method
 - Dual line method
 - Triple line method
 - Single ellipse method
 - Single ellipse and single line method, Trace&H
- Angle
 - PolyLine,TwoLine
- %Stenosis
 - Diameter method
 - Square meter method
- A/B ratio
 - Diameter ratio
 - Square meter ratio

8.1.2 Generic Measurement in CM mode

- CFVP
 - point
 - profile

8.1.3 Generic Measurement in M mode

- Depth
- Distance

- Time
- Slop (Velocity)
- Heart rate
- %Stenosis
- A/B ratio
 - Diameter ratio
 - Time ratio
 - Speed ratio

8.1.4 Measurement in PW mode

- Velocity
- Time
- Acceleration
- PS (Peak Speed in systole period)
- ED (The speed in the end of diastole period)
- MD (Minimum speed in diastole period)
- Mean Vel(Max)
- Mean Vel(Mean)
- PI (Pulsatility Index)
- RI (Resistance Index)
- PS and ED ratio
- ED and PS ratio
- A and B ratio (A/B ratio)
 - Speed ratio
 - Time ratio
- Auto Flow Volume(optional)
- MaxPG (maximum pressure gradient)
 - MeanPG (Mean pressure gradient)
 - SV (Stroke Volume)
 - Each volume diameter cardiac
 - Time mean speed in each

stroke volume

- Cardiac output
- Heart rate
- $SV(LVOT)/SV(RVOT)$

8.2 Measurement in different applications

8.2.1 Abdominal Measurement

- General abdomen
- Difficult abdomen
- Kidney
- Renal vessel
- Abdominal trauma

8.2.2 Small Part Measurement

- Thyroid
- Breast
- Testis
- Musculoskeletal
- Upper and lower extremity joint
- Nerve block

8.2.3 Vessel Measurement

- Carotid artery
- Upper artery
- Upper vein
- Lower artery
- Lower vein
- Vessel puncture
- Transcranial Doppler

8.2.4 Gynecology Measurement

- Uterus and Pelvic
- Follicle

8.2.5 Urology Measurement

- Bladder
- Prostate
- Renal
- Kidney and ureter
- Pelvic Floor dysfunction

8.2.6 Pediatric Measurement

- Neonatal Head
- Neonatal Abdomen
- Pediatric Abdomen
- Pediatric Hip

8.2.7 Obstetrics Measurement

- OB Early
- OB Mid
- OB Late
- Fetal Heart

8.2.8 Cardiac Measurement

- General
- LV
- MV
- Ao
- LA
- RV
- TV
- PA

- RA

8.2.9 Auto NT (Nuchal Translucency) measurement(optional)

- Automatically detect Nuchal Translucency in interest box
- Automatically report thickness result of NT

8.2.10 Auto IMT (Intima-Media Thickness) measurement

- Automatically detect intima media thickness in interest box
- Automatically report the result of IMT
- Available in linear probe

8.2.11 Live IMT (Intima-Media Thickness) measurement (optional)

- Real-time automatically display IMT items with the different ROI positions
- The IMT items include: max, min, average, SD, points (how many points are used for the result), size of ROI .
- Available in carotid application

8.2.12 Auto IT (Intracranial translucency) measurement(optional)

- Support Auto IT(Intracranial translucency) measurement
- Draw the ROI and the system analyses and displays the result

8.2.13 Auto Follicle(2D/3D)(3D optional)

- Just click on the area of follicle in B mode, the area of this follicle will be reported automatically
- Report the area of different follicle in the volume data automatically

8.2.14 Smart 3D Volume Measurement(optional)

- Trace the margin of the irregular circle in different slices of volume data in irregular shape
- Automatically report the volume of the irregular object

8.2.15 VAim OB

measurement(optional)

- VAim OB is an intelligent tool for fetal growth calculation, just one touch to activate the measurement items (BPD, OFD, HC, AC, FL, HL) and get the results, helps to make clinical decisions quickly and confidently, improving the speed and ease of exams
- The intelligent results will be add into the worksheet and report automatically

8.2.16 VAim Hip measurement(optional)

VAim Hip is an intelligent solution in the assessment of DDH(Developmental Dysplasia of Hip) with one simple touch.

- Based on ‘Ped HIP’ application

8.2.17 VAim Follicle(2D)measurement (optional)

An intelligent tool for follicle calculation, one touch to get the follicle status, dedicated for women’s reproductive health.

- Choose left or right follicle
- Automatically identity all the follicles with different colors and calculate follicle volume and diameter

8.2.18 VAid Breast(optional)

VAid Breast is an automatic tool for breast lesion detection in real-time or on stored images (static & cine)

- For static image: Depicts boundaries of the Breast lesions and displays the size
- For cine(real time scan or stored cine):The number and position of the Breast lesion can be indicated in real time.

8.2.19 VAid Thyroid(optional)

VAid Thyroid is a tool for Thyroid lesion detection in real-time or on stored images (static & cine)

- For static image: Depicts boundaries of the thyroid lesions and displays the size
- For cine(real time scan or stored cine):The number and position of the Thyroid lesion can be indicated in real time.