The BDC system can be easily installed in any mobile X-ray Chest Survey Cabinet, mounted on truck.

Standard lay-out of the screening room, to optimize the flow of patients.

Patients room (1) separated from the personnel working area, to protect technicians from contagious diseases.


Working area equipped with working station (2), wardrobe, sofa, air conditioning system (3), small table (4).

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Data subjected to alteration without prior notice. DBC series - 02/2003

BMI

ESSE 3 Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.4ci.com

BMI

BCD Digital Chest Screening Radiographic System, dedicated to mass chest examinations for early diagnosis of tuberculosis.

High-quality digital chest images with very low radiation dose delivered to the patient (dose level below 5microSv).


Chest images realized by the horizontal movement of a sharp X-ray beam synchronized with the movements of the X-ray detectors. Complete chest scan in less than 4 seconds.

Linear X-ray detectors with extremely high contrast sensitivity due to the special X-ray beam

collimation, reducing the influence of the harmful scattered radiation.

High spatial resolution digital images with a very small focal spot X-ray tube.

Detectors array composed of 1280 photodiodes covered by CsI scintillating screen Xray detector.

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
BCD Digital Chest Screening System is composed of:

- X-ray protected cabin, with motorized footrest to adjust patient position.
- H.F. X-ray Generator with rotating anode X-ray tube.
- The two **BDC** working stations consists of the technician's and the radiologist's computer.
- **BDC** System uses **DICOM 3.0** software compatible to and easy to integrate to a RIS.

The 3 mm lead equivalent protection of the cabin walls and windows makes it possible to install the **BDC** system anywhere with no need of special additional X-ray protection.

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