



BESDATA

Video Laryngoscope





- * **One-button operation**
- * **Automatic power off**
- * **Smart Anti-fog technology**
- * **250 minutes battery with screen low battery alarming notice**
- * **AV output to connect to monitor for training usage**



Technical Data

Video laryngoscope handle: Size 180x70x110mm
Weight 180gr
Angle of vision $\geq 60^\circ$

Battery: Power 3.7V lithium battery
Discharging time of battery >1hr
Charging time >300 times
Resolution $\geq 3.72\text{lp/mm}$

Monitor: Light source high-intensity LED
LCD Display : **BESDATA 3** : 3" and **BESDATA 7** : 7" inches color display
Camera CMOS
Rotating angles of the monitor in front to back $0^\circ \sim 130^\circ$
Rotating angles of the monitor in left and right $0^\circ \sim 270^\circ$

Video laryngoscope reusable handles

Material: Anti-Fog medical grade polycarbonate
Blades size: MAC-Style 3 size for Adult/Pediatric/Neonate



Video laryngoscope disposable blades

Material: Anti-Fog medical grade polycarbonate
Blades size: MAC-Style 1 size for Adult
Packaging: Packaged sterile for single use



Video laryngoscope metal blades

Material: 316 stainless steel
Blades size: D-style 1 size adult for difficult airway



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

