



Comprehensive and Powerful Video GASTRO & COLON VIDEOSCOPIES




ES-EGS69 / ES-EGS70
Gastrointestinal & Colono VideoScope

 **ESSE3** srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail: danielec@esse3-medical.com
web: www.esse3.dreamgest.com



Gastrointestinal & Colono VideoScope



Item/Model	Gastrointestinal videoscope ES-EGS69 / ES-EGS691	Colonovideoscope ES-EGS70 / ES-EGS701
Field of view	120° Forward viewing	120° Forward viewing
Depth of View Field	5mm~100mm Fixed Focus	5mm~100mm Fixed Focus
Range of Tip Bending	Up 210° , down 100° Left 100° , right 100°	Up 180° , down 150° Left 150° , right 150°
Outer Diameter of Distal Fnc.	ø 9.8mm/ø 9mm	ø 13.2mm/ø 14.3mm
Outer Diameter of Insertion Tube	ø 9.8mm/ø 9mm	ø 12.8mm/ø 13.7mm
Inner Diameter of Instrument Channel	ø 2.8mm/ø 2.2mm	ø 3.8mm/ø 4.2mm
Working Length	1025mm/1025mm	1300mm/1500mm/1680mm
Waterproof	•/•	•/•

Standard Configurations:

- Gastrointestinal and Colono VideoScope main unit*1;
- Leak detection cover*1;
- Plug the mouth cap*5;
- Hand pressure leak detector*1;
- Cleaning brush*1;
- Biopsy forceps*1;
- Water bottle assembly*1;
- Mouth pad*2;
- Technical Instructions*1;
- User Manual*1;
- Carton box*1



*Specification design and accessories are subject to change without any notice on obligation on the part of manufacturer.