

Endoscope Storage and disinfection Cabinet



Features

- UV Wavelength: $\lambda = 253.7\text{nm}$
- UV Radiance: $> 90\text{uw/cm}^2$
- UV Radiation: $300,000\text{uw.s/cm}^2$
- UV lamp lifetime : 7500hours
- Disinfection volume: 99.9%
- Power supply : 220V 50Hz (110V, 60Hz as optional)
- Power: $< 200\text{Watts}$

With integrated air moisture removing system, UV radiation system makes the carbine realizes the best storage and disinfection surrounding for the flexible endoscope.

With integrated microcomputer control, the air disinfection and moisture time can be adjusted.

With the time-set system, the nurses can set the time to disinfect and circle ending time, we recommend 20 minutes a circle.

With the hook for the flexible forceps , maximum to 30 pieces forceps can be stored.

Specification

Item-No.	Description	Size	Endoscope storage
ESDC-100	Single door	2000x680x530mm (HxWxD)	4 pieces
ESDC-200	Double door	2000x1200x530mm (HxWxD)	8 pieces