





MODEL: **ES-EH58**



Left 170° Right 170°
Screen rotation angle

MODEL	ES-EH58	ES-EH28	ES-ES28
Field of View	120°	120°	120°
View Direction	0°(Direct-Vision)	0°(Direct-Vision)	0°(Direct-Vision)
Depth of View	2-50mm	2-50mm	2-50mm
Working Length	600mm	700mm	600mm
Insertion Out Dia	15Fr(Φ5mm)	8.5Fr(Φ2.85mm)	8.5Fr(Φ2.85mm)
Biopsy Channel	6Fr(Φ2.0mm)	3.6Fr(Φ1.2mm)	3.6Fr(Φ1.2mm)
Deflection	Up: 170° Down: 170°	Up: 275° Down: 275°	Up: 180° Down: 130°



Applications:

Model: **ES-EH58**

ES-EH58 is used in human medical as a Nasopharyngoscope.

Model: **ES-EH28**

ES-EH28 is a Ureteroscope used for urinary stones and kidney stones in the urology department.

Model: **ES-ES28**

ES-ES28 is a pediatric bronchoscope used in the respiratory department.



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail: danielec@esse3-medical.com
web: www.esse3.dreamgest.com

