



ESMC-340

Medical Air Compressor





ESMC-340

Features

- CE marked, stable performance, safe and reliable
- Compact design and easy to transfer
- Assembled with clean control module, supplying cleaner and safer air to the patients
- High effective sound-absorbing materials, noise reduction
- Equipped with various ventilator of famous brand in the world

Technical Data

- Output pressure: 0.4Mpa
- Output flow (continuous): $\geq 40\text{L/min}$
- Output flow (peak): $\geq 75\text{L/min}$
- Peak flow: $\geq 180\text{LPM}$
- Noise(Actually): $< 52\text{dB(A)}$



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.esse3-medical.com

