



DATA SHEET 

Model: ES-PL

Mechanical plasma extractor



- Manual system – accept all kinds of blood bags.
- Frame and construction in stainless steel
- Transparent plate for visual control red cells / plasma
- Powerful spring
- Dimensions (WxDxH) : 17x26x22 cm
- Gross weight: 3,2 kg



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com



CE