



ESCH-105L2

Low speed inner water channel handpiece



Technical specification

Working pressure	0.3mpa – 0.35mpa
Applicable burs	2.335 – 2.355mm
Rotation speed	20000 – 30000rpm
Resistant to high temperature Disinfection	135°C
Noise sound	≤ 70dB



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail: danielec@esse3-medical.com
web: www.esse3-medical.com

