



Cod: **INTRA SURGERY 2**

IMPLANTOLOGY MICROMOTOR TRAUSS

Speed: 20:1 from 30 to 2000 rpm / 32:1 from 20 to 1250 rpm Torque:  
20:1 from 5,00 to 55 Ncm / 32:1 from 5,00 to 65 Ncm Gear ration: 1:5 /  
1:4 / 1:1 / 16:1 / 20:1 / 27:1 / 32:1 / 64:1


Handpiece: Traus CRB26LX (with optical fiber) Traus CRB26XX (without optical fiber)

Torque maximum: 80 Ncm

Motor e-type: Traus MBP10SL (sx) standard connection

Pump maximum flow: 75 ml/min.

- 10 programs with memory function
- Protection of overload
- Torque and speed shown on the display when the handpiece is working so as to control the phases of the operation
- Foot Control Ergonomic



CE



- |                         |                              |
|-------------------------|------------------------------|
| 1 - Display volume      | 4 - Display direction        |
| 2 - Display program     | 5 - Display speed/torque     |
| 3 - Display gear ration | 6 - Display optical function |

LED handpiece 20:1 Traus - CRB26LX



LED handpiece 32:1 Traus - CRB27LX



## HOSE SET FOR IMPLANTOLOGY MOTOR



CE

item: **575/43**

Spare throwaway sterile tubes for IntraSurgery 2.0.  
Pack 10 pcs

item: **554/43**

Spare throwaway sterile tubes for surgical units.  
Pack 10 pcs

item: **565/43**

Spare throwaway sterile tubes for UltraSurgery.  
Pack 10 pcs



CE

Cod: **INTRA SURGERY 2**

MICROMOTEUR D'IMPLANTOLOGIE TRAUSS

Vitesse: 20:1 de 30 à 2000 rpm / 32:1 de 20 à 1250 rpm

Couple: 20:1 de 5,00 à 55 Ncm / 32:1 de 5,00 à 65 Ncm

Ration des engrenages : 1:5 / 1:4 / 1:1 / 16:1 / 20:1 / 27:1 / 32:1 / 64:1

Pièce à main: Traus CRB26LX (avec fibres optiques) Traus CRB26XX (sans fibres optiques)

Torque maximum: 80 Ncm

Moteur e-type: Traus MBP10SL (sx) connection standard

Flux maximal de la pompe: 75 ml/min.

- 10 programmes avec fonction de mémoire
- Protection de surcharge
- Couple et vitesse affichées sur l'écran quand la pièce à main est en fonction pour contrôler les phases de l'opération
- Contrôle au pied ergonomique



Pièce à main LED 20:1 Traus - CRB26LX



Pièce à main LED 32:1 Traus - CRB27LX



- |                                    |                            |
|------------------------------------|----------------------------|
| 1 - Ecran volume                   | 4 - Ecran de direction     |
| 2 - Ecran programme                | 5 - Ecran vitesse/couple   |
| 3 - Ecran de ration des engrenages | 6 - Ecran fonction optique |

## JEU DE TUYAUX POUR LE MOTEUR D'IMPLANTO@CGIE



CE

item: **575/43**

Tuyaux de rechange stériles jetables pour IntraSurgery 2.0.

Pack 10 pcs

item: **554/43**

Tuyaux de rechange stériles jetables pour unités chirurgical.

Pack 10 pcs

item: **565/43**

Tuyaux de réchange stériles jetables pour UltraSurgery.

Pack 10 pcs