




ESCH-3000

PORTABLE DENTAL COMPRESSOR WITH DELIVERY SYSTEM

"The pictures are for reference only. The configuration is subject to the actual product."



Instrument classification: Class I

Power Source: Electricity

Material: metal, Plastic

Quality Certification: CCC

Portable Dental Turbine Unit with Air Compressor Free

Handpiece tubing: 2pcs

Warranty: One Year

Three Way Syringe: Yes

Saliva Ejector: Yes

Foot Pedal: Yes

Feature: with Air Compressor

Package: Carton

Air bottle: 13L

Net weight: 48 kg

Gross Weight: 65 kg

Single package size: 57X75X107 cm



Dental Air Pump

Power supply: 110V/220V, AC 60Hz/50Hz

Current: 1.6A

Power: 550W

Gross weight: 9.5kg/13Kg

Package size: 29*18*24cm

Package size: 40*32*27cm



Options Accessories

- 1.High-speed Dental handpiece(2-hole/4-hole)
- 2.Low-speed Dental handpiece(2-hole/4-hole)
- 3.Built-in Electric Motor
- 4.LED Sensor Lamp
- 5.Endoscope System Wifi Mode
- 6.Built-in Fiber Optic System(6-hole)



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail: danielec@esse3-medical.com
web: www.esse3-medical.com

