



# ESX-60P

Dental x-ray  
portable unit



DENTISTRY



# ESX-60P

Dental x-ray portable unit



Cod. **ESX-60P**



- This portable dental X-ray unit is with high frequency.
- Small body, light weight and nearly no radiation.
- LCD display and excellent image quality.
- Portable carrying case, save more space.
- Touch buttons on the control panel
- Rechargeable batteries and charger included.
- Suitable unit mainly for oral pre-treatment for learning internal organizational structure, root depth and so on, is indispensable in daily life clinic equipment, especially for dental implant surgery.
- It can be also connected with a sensor.

• **Technical features:**

Tube voltage	60kV
Tube current	1mA
Exposure time	0.02~2
Frequency	30kHz
Tube focus	0,3 x 0,3
Distance from skin to cone	130 mm
Rechargeable battery	DC14.8V, 6400mA
Rated power	60VA
Charger input voltage	AC100V-240V ±10%
Charger Output voltage	DC 16,8 V
Ambient temperature	10-40 °C
Relative humidity	10%-75% RH
Net weight	2,5kg
Size	138x165x185 mm
Packing size	220x240x275 mm



Control panel



Carrying case



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

