

ES-BSEAL



VESSEL SEALING SYSTEM




Sealer Modes :

Mode	Max. Power Watt.	Max. Short Circuit Current	Peak to Peak Voltage	Operating Frequency	Crest Factor
Soft	150	4.4 Amps	560 (280)	480 khz	1.5
Force	150	4.4 Amps	700 (350)	480 khz	1.5

Bipolar Modes :

Mode	Max. Power Watt.	Max. Short Circuit Current	Peak to Peak Voltage	Operating Frequency	Crest Factor
Macro	120	2.2 Amps	560 (280)	480 khz	1.5
Macro	80	2.2 Amps	560 (280)	480 khz	1.5
Bi-Cut	120	2.2 Amps	620 (310)	480 khz	1.5

Display	: 7 segment LED
Dimension (W x H x D)	: 290mm x 105mm x 260mm
Weight	: 4 kg
Keyboard	: Feather Touch
Manufactured and Tested acc. to std.	: IEC 601-1 / EN 60601-1
Input Supply Voltage	: 230 VAC + 15% @ 50/60 Hz (110 V optional)



Standard Accessories

- Single Pedal Foot Switch
- Double Pedal Foot Switch
- Autoclavable Bipolar Forcep with Cord (At Additional Price)
- Laparoscopic Sealer Instrument (At Additional Price)



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.4ci.com

