





# ES-MFM-CNS

Central Monitoring System



INTENSIVE CARES



ES-MFM-CNS is a comprehensive central nursing system for the obstetric cares. It is an integrated system that covers the entire continuum of obstetric care across the whole pregnancy stage from the first antepartum visit to post delivery care. The system displays data and information transmitted from fetal and maternal monitors and generates alerts for critical events. Easy access to well-organized information improves efficiency, allowing you to focus more attention on patients and achieve the best quality of care.

- Connection of up to 128 bedside fetal/maternal monitors
- Bi-directional communication of patient ID and alarming setup
- Time synchronization
- Built-in antepartum CTG & NICHD analysis function
- Up to 200 groups of NIBP data for each record
- Remote monitoring from different locations via CNS Client
- Supports for 32 people login at the same time
- Built-in partogram for labor stage management
- Export PDF&JPG format report
- Windows 10 Compatibility



ESSE3 srl, Via Garibaldi 30  
 14022 Castelnuovo D.B. (AT)  
 Tel +39 011 99 27 706  
 Fax +39 011 99 27 506  
 e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
 web: [www.esse3.dreamgest.com](http://www.esse3.dreamgest.com)

