



ESTX-D3

Benchtop cytocentrifuge



LABORATORY



ESHL!8'

.....6YbW hcd W h c W b h f] Z i [Y



Cod. ESTX-D3

Benchtop centrifuge for cytological specimens that transfer cells from fluids to glass slides for microscopic examination.

Microprocessor control and powered lid locking.

Control panel significantly facilitates daily laboratory routine. Digital display. With 25 programs.

TECHNICAL FEATURES	
Cod.	ESTX-D3
N. of programs	25
Max speed	3000 r/min
Max RCF	730×g
Volume of single tube	4 ml
Qty of total tubes	6
Programmable Timer	1 ~ 99 min
Motor	DC brushless motor, Microprocessor control
Power supply	AC220V 50Hz 5A
Dimensions	370×432×250 mm (L×W×H)
Weight	14 kg



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

