



# ES-MICROPHOR

Electrophoresis reader



LABORATORY

# ES-MICROPHOR

Electrophoresis reader



## Technical features

- ◆ LED light source and high gain fototransistor.
- ◆ Reading system in absorbance with self-regulation of the source bright.
- ◆ Reading media in cellulose acetate (trasparent and opaque)
- ◆ Reading media in agarose gel
- ◆ 8 reading samples at a time.
- ◆ Connection to PC via uSB
- ◆ External power supply 24v.
- ◆ Dimension: 11.6x15.6x6.3 cm (LxPxH)



## Software features

- ◆ Pre-programmed Analysis: Serum proteins, Heamoglobins, Lipoproteins.
- ◆ Possibility to define adicional exams.
- ◆ Automatic processing of electrophoretic pattern.
- ◆ Possibility of correction of paths
- ◆ Personal management and printing reports.
- ◆ Report header customization
- ◆ Operating system: Microsoft windows compatible.



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706

esse3@chierinet.it  
www.esse3-medical.com

