



EG.66G!J % \$\$\$ / EG.66G!J % \$\$\$` Vertical laminar flow hoods



LABORATORY



ES!66G!J % \$\$\$#% \$\$

Vertical laminar flow hoods




LCD Display

Airflow speed, UV timer, UV work time, system work time, real time



Waterproof Socket

2 waterproof sockets in the side panel, for to use small devices inside the cabinet



UV lamp

Emission of 253,7nm for an efficient decontamination




Gas Tap



Side Glass Window

The transparent side glass windows maximize light and visibility inside the cabinet,

providing a bright and open working environment.



Technical features:

Model	ES-BBS-V1300	ES-BBS-V1800
External Size (WxDxH)	1310x750x2000mm	1800x700x2010mm
Internal Size (WxDxH)	1200x645x610mm	1700x645x610mm
Working Height	720mm	720mm
Display	LCD Display, with air flow display and real time	LCD Display, with air flow display and real time
Airflow Speed	Average of 0,3~0,5m/s	Average of 0,3~0,5m/s
Material	Main Body: Cold-rolled steel with anti-bacteria powder coating Work Table: 304 stainless steel	Main Body: Cold-rolled steel with anti-bacteria powder coating Work Table: 304 stainless steel
Pre-filter	Polyester fiber, washable	Polyester fiber, washable
HEPA Filter	99,999% efficiency at 0,3um	99,999% efficiency at 0,3um
Noise	<60dB	<60dB
Front Window	Motorized, Side and Front Windows: 5mm toughened glass, anti UV	Motorized, Side and Front Windows: 5mm toughened glass, anti UV
Max Opening	520mm	520mm
Fluorescent Lamp	28Wx1	36Wx1
UV Lamp	Emission of 253,7 nanometers	Emission of 253,7 nanometers
Consumption	400W	450W
Waterproof socket	Two, total load ≤ 500W	Two, total load ≤ 500W
caster	Universal casters with leveling feet	Universal casters with leveling feet
Power Supply	110/220V±10%; 60/50Hz	110/220V±10%; 60/50Hz
Standard Accessories	Fluorescent lamp, UV lamp x2, Base stand, Tap, Waterproof socket x2	Fluorescent lamp, UV lamp x2, Base stand, Tap, Waterproof socket x2
Gross Weight	170kg	255kg
Packing Size (WxDxH)	1470x980x1570mm	1870x980x1570mm



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

