



OL135B Series

Binocular Microscope



MICROSCOPY



OL135B Series

Binocular microscope

Description

Model OL135B / OL135BL400 / OL135BL:

- Binocular head 30 ° inclined, 360 ° rotatable.
- Diopter adjustment ± 5 (diopters) and interpupillary distance.
- Eyepiece WF10x / 18mm.
- DIN 160mm 4x / 0.10 - 10x / 0.20 - 40x / 0.65 - 100x / 1.25 Oil objectives. 100x not supplied for model OL135BL400
- Nosepiece nosepiece 4 places (click stop).
- Revolver tilted towards the stand to facilitate the insertion of the slides
- Abbe condenser A.N. 1.25 iris diaphragm / filter holder height adjustable by knob.
- 125x115mm double-layer translating table with coaxial controls.
- Coaxial coarse / micrometric focusing with clutch and limit switch.
- Focus intensity adjustment system and table height limit switch system to avoid accidental breakage with the glass sample holder slide.
- Adjustable 12V / 20W halogen lighting (OL135B).
- 0.6 W LED lighting adjustable with rechargeable batteries (OL135BL).



Lighting

Choice of Halogen or LED.
Halogen lighting 12V 20W
LED lighting 0.6



Focus

Coaxial macro-micrometric focusing with limit switch and friction adjustment system. Double-Layer type translating table with coaxial controls on the left



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

