



ESPX-B Series

Biochemistry Incubators



LABORATORY



Features

- * ESPX-BI series: LED display.
- * ESPX-BII & BIII series: LCD display.
- * Built-in printer for ESPX-BIII series.
- * Polished stainless steel chamber.
- * Round angle structure , easy to clean.
- * Universal Castors for easy movement.
- * CFC-free eco-friendly refrigerant.
- * Microprocessor PID temperature control, accurate and reliable.
- * Power socket integrated inside for additional small equipment.

Optional accessories:

- * Test port

Safety function

- * Temperature protection threshold value can be set.
- * Auto-off protection and alarm for abnormal situation.
- * Parameters can be saved automatically when power off or system error.

Technical parameters

Model	ESPX-B80I ESPX-B80II	ESPX-B150I ESPX-B150II III	ESPX-B200I ESPX-B200II III	ESPX-B250I ESPX-B250II III	ESPX-B300I ESPX-B300II III	ESPX-B400I ESPX-B400II III
Capacity	80L	160L	202L	251L	304L	397L
Temp. Range	5~50 °C(type II 0~65°C)					
Temp. Fluctuation	±1°C (type II & III: heating: ±0.5°C to 1°C);					
Temp. Uniformity	±2 °C					
Shelves	2pcs/adjustable	3 pcs/adjustable				
Time Range	0~9999min					
Refrigerant	R134a (CFC free, eco-friendly)					
BOD socket	1pc					
Consumption	650W	710W	760W	810W	960W	1280W
Power Supply	AC220V±10%, 50/60Hz					
Internal Size(mm)	400*370*550	450*420*850	450*450*1000	480*490*1070	520*500*1170	580*540*1270
External Size(mm)	540*570*1070	600*670*1470	590*640*1520	620*680*1590	660*690*1690	720*740*1790
Package Size(mm)	700*730*1260	740*810*1650	750*800*1710	780*840*1780	820*850*1880	880*900*1980
Gross Weight(kg)	85	145	146	147	157	170



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

