



ES-709/O



- 1 Outlet supply / Alimentation de sortie
- 2 Overload fuses / Fusibles de protection
- 3 Main switch / Interrupteur général



HORIZONTAL SHAKING STIRRER

The horizontal shaking stirrer **ES-709/O** has a work surface with horizontal moving of 30 mm. It is suitable for holding and fixing, by adjustable clips, beakers and flasks of different sizes. Electronic regulation of speed from 50 to 250 oscillation per minute. Possibility of continuous or pre-set working from 1 to 999 minutes. Possibility to memorize 4 disjoint programming channels of time and speed values.

Polycarbonate anti-scratch panel with visualization of speed (RPM) and remaining time to finish cycle.

Outer casing in powdered painted steel acid resistant and work surface in stainless steel AISI 304 2B glazed with coverage of not slip material. Main switch bright.

TECHNICAL SPECIFICATIONS

Overall dimensions WxDxH: 630x455x310 mm.

Work surface dimensions WxD: 550x 330 mm.

Weight: 31 Kg.

Max capacity of work surface: 15 Kg.

ELECTRICAL SPECIFICATIONS

Power supply: 230 V - 50 Hz.

Electrical input: 70 W

Overload fuses: 2 fuses x 0.5 AF (5 x 20) mm.



AGITATEUR OSCILLANT

HORIZONTAL

L'agitateur oscillant horizontal **ES-709/O** a une surface de travail avec mouvement horizontal de 30 mm. Il est indiqué pour l'agitation, par ses clips réglables, de béciers et de flacons de tailles différents. Régulation électronique de la vitesse de 50 à 250 oscillations par minute.

Possibilité de fonctionnement continu ou pré-réglé de 1 à 999 minutes. Possibilité de mémoriser 4 canaux de programmation disjointes de valeurs de temps et vitesse.

Panneau anti-rayures en polycarbonate avec visualisation de la vitesse (TPM) et du temps restant pour terminer le cycle.

Boîtier externe en acier peint époxy résistant aux acides et surface de travail en acier inox AISI 304 2B vitré avec une couverture de matériel non glissant. Interrupteur général lumineux

SPECIFICATIONS TECHNIQUES

Dimensions externes LxPxH: 630x455x310 mm.

Dimensions surface de travail LxP : 550x330 mm.

Poids : 31 kg

Capacité maximale de la surface de travail : 15 kg.

SPECIFICATIONS ELECTRIQUES

Alimentation : 230 V – 50 Hz.

Puissance électrique : 40W

Fusibles de protection : 2 x 0.5 A F (5x20) mm.



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)

tel +39 011 99 27 706

fax +39 011 99 27 806

e-mail esse3@chierinet.it

web : www.esse3-medical.com

