



# ES-JOYM

Electronic baby scale  
with height rod



LITTLE ACCESSORIES



# ES-JOYM

Electronic baby scale with height rod

Cod. **ES-JOYM**



| Technical Features          |   |
|-----------------------------|---|
| <b>Max. capacity</b>        | 20kg  |
| <b>Division</b>             | 5/10g   |
| <b>Structure</b>            | Plastic ABS structure with washable tray  |
| <b>Functions</b>            | Tare, Hold weight, Unit weight, Weight memory and difference calculation, External calibration, Programmable auto-OFF |
| <b>Display</b>              | Large LCD 24mm 4 digit, easy to read  |
| <b>Power</b>                | 1 alkaline battery (PP3 type) or external power adapter   |
| <b>Weight</b>               | 3 kg  |
| <b>Dimensions (WxLxH)</b>   | 305 x 520 x 120 mm  |
| <b>Included Accessories</b> | Power adapter, Low batteries indicator  |
| <b>Optional Accessories</b> | Alkaline battery (PP3 type)   |

  

| Cod.           | Descr.  |
|----------------|---|
| <b>ES-JOYM</b> | Baby scale with mechanic height rod on the tray |
| <b>ES-JOYN</b> | Baby scale with digital height rod on the tray  |



ESSE3 srl, Via Garibaldi 30  
 14022 Castelnuovo D.B. (AT)  
 Tel +39 011 99 27 706  
 Fax +39 011 99 27 506  
 e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
 web: [www.esse3-medical.com](http://www.esse3-medical.com)

