



ESF10 Line BABY SCALES



LITTLE ACCESSORIES



ESF10 Line BABY SCALES

Cod. **ESF10-7550** - ELECTRONIC BABY SCALE



Directive 93/42/CEE

Max. capacity	20 Kg
Min. capacity	30 g
Accuracy	10 g
Weight unit	Kg / lb
Overload	0 Ld
Battery level indicator	Yes

Cod. **ESF10-7800** - ELECTRONIC BABY SCALE



Max. capacity	20 Kg
Accuracy	5 g
Weight unit	Kg / lb
Overload	0 Ld
Battery level indicator	Yes

Cod. **ESF10-7660** - ELECTRONIC BABY SCALE



Max. capacity	20 Kg
Accuracy	5 g
Weight unit	Kg / lb
Overload	0 Ld
Battery level indicator	Yes

ESF10 Line BABY SCALES

Cod. **ESF10-725** - MECHANICAL BABY SCALES WITH SLIDING WEIGHTS

Tray dimensions : 57 x 30 x 9 cm

Max. capacity	20 Kg
Accuracy	10 g



Cod. **ESF10-7250** - MECHANICAL BABY SCALES WITH SLIDING WEIGHTS

Max. capacity	16 Kg
Accuracy	10 g



Cod. **ESF10-7463** - MECHANICAL BABY SCALES WITH SLIDING WEIGHTS

Stainless steel tray dimensions : 52 x 29 x 5 cm
Accuracy class : III

Max. capacity	16 Kg
Accuracy	10 g





ESF10 Line BABY SCALES

Cod. **ESF10-7465** - ANALOGIC BABY SCALE



Stainless steel tray dimensions : 52 x 30 x 5 cm

Max. capacity	20 Kg
Accuracy	50 g / 10 kg 100 g / 20 kg

Cod. **ESF10-7600** - ANALOGIC BABY SCALE



Stainless steel tray dimensions : 52 x 30 x 5 cm
Accuracy class : IIII

Max. capacity	20 Kg
Accuracy	10 kg / 50 g 20 kg / 100 g

Cod. **ESF10-501** - SALTER BABY SCALE WITH LARGE DIAL



With baby trouser

Max. capacity	25 Kg
Accuracy	100 g
Optional accessories	
ESF10-501/10	Baby trouser
ESF10-501/20	Carrying case

