



# OTOSCOPE

**Fiber optic halogen light.**

***Lumière halogène à fibres optiques.***



Cod. **ESF8-660019**



- 2,5 V Fiber optic halogen light.
- Otoscope head: pin contact fitting.
- Fixed lens 4X in plastic mount may be removed for introduction of instruments.
- Plastic tube with 7 ear disposable specula.
- Battery handle: made in plastic material for 2 AA batteries.
- Replacement bulb: 2,5 V.

Optional Accessories:

Disposable specula - 2,5 or 4 mm



- *Lumière halogène à fibres optiques de 2,5V*
- *Tête d'otoscope: fermeture à vis.*
- *Le lentille 4X sur monture plastique peut être enlevé.*
- *Tube en plastique et 7 spéculums à usage unique.*
- *Manche à piles: Fait en plastique pour 2 piles type AA.*
- *Ampoule de rechange: 2,5 V.*

*Accessoires Optionnels:*

*Spéculums à usage unique - 2,5 ou 4 mm*



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

