

### ESKD-202A-7 Head light



- ◆ With Built-in Battery .
- ◆ Two battery , one is 4h, another one is 2h.
- ◆ Feature: 3W LED illumination, sunlight color, cold light, cableless, two batteries
- ◆ Function: Quick of battery replacement, brightness adjustable, can be used with magnifiers
- ◆ Variable function: Variable color on the surface
- ◆ Application: Surgery room, stomatology, ENT


Technical Data									
Voltage	Power	C.T.	Illumination	Size of light spot	DC Time	Lamp Life	Charge Time	Weight	Package
AC 95-245V	LED 3W	5500K	60000LX	10-110mm	4h/2h	50000h	5h	360g	205*190 *127mm

### ESKD-202A Headlight



- ◆ With Built-in Battery .
- ◆ Two battery , both 4h
- ◆ Brightness adjustable
- ◆ Adapt top world LED light source, color rendering index CRI greater than 90%
- ◆ Illumination spot is clear for excellent optical design
- ◆ Cableless, two batteries, the battery energy capacity is not attenuated after the battery repeated be charged and discharge 500 times, quick of battery replacement, brightness adjustable, can be used with magnifiers
- ◆ Application: Cardiothoracic surgery, neurosurgery, ENT, general surgery, gynecology, dentistry

Technical Data									
Voltage	Power	C.T.	Illumination	Size of light spot	DC Time	Lamp Life	Charge Time	Weight	Package
AC 95-245V	LED 5W	5500K	70000LX	10-110mm	4h/4h	50000h	5h	420g	270*250 *160mm



### ESKD-202-1 Head loupe



LED light optional



#### Technical Data:

1. Voltage: AC95-245V
2. Lamp power: LED 3W
3. C.T.: 5500±500K
4. Lamp life: 50000h
5. Illumination: 65000lx
6. Size of light spot: 60mm
7. N.W.: 0.26KG
8. DC Charging: 5h
9. DC Power supply: 4h
10. Package size: 170\*125\*52mm

- ◆ Feature: 3W LED illumination, sunlight color, cold light, Simple, Economical, convenient to carry
- ◆ Function: Used with magnifier
- ◆ Application: Outpatient service, stomatology, ENT

### ESFD-503G Binocular Loupe



#### Characteristics:

There are 2.5x, 3.5x 4.0x 5x 6x 6.5x

Galileo magnifying glass;	2.5X	3.5x	4.x	5x	6x	6.5x
working distance:	420mm;	420mm	420mm	420mm	420mm	420mm
field of view:	130mm;	100mm	80mm	70mm	60mm	55mm

glasses frame type; spiral adjustment of interpupillary distance

### ES-HLD Binocular Loupe



There are 2.5x, 3.5x 4.0x 5x 6x 6.5x

Galileo magnifying glass;	2.5X	3.5x	4.x	5x	6x	6.5x
working distance:	420mm;	420mm	420mm	420mm	420mm	420mm
field of view:	130mm;	100mm	80mm	70mm	60mm	55mm

glasses frame type; spiral adjustment of interpupillary distance