



EMERGENCY BAG



Cod. 9G: -!, - \$&)

Dimensions:
Width 47 cm
Depth 26 cm
Height 30 cm

Compact and capable emergency bag made of strong waterproof red fabric with zipper closure on top and non-slip feet for stability, with reflective band, transportable with solid handles or sling with a snap hook.

In the inner part, the provided dividers allow to organize the space to storage as needed for emergency response.

The bag is equipped of containers with zipper closure at the four corners. The container in the frontal part is equipped of additional housings useful for storing equipment.

Easy cleaning.



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

