



ESF9 Line

Emergency equipment



LITTLE ACCESSORIES



ESF9 Line

Emergency Equipment



Silicon Reanimation Bag

100% latex free. Completely autoclavable. Made of pure, durable, silicon.
 All parts can be easily removed to ease cleaning or disinfecting.
 Complete with pressure release valve to minimize the risk of over inflating the lungs.
 The release valve can be locked if needed.



ESF9-0720S



ESF9-0721S



ESF9-0722S

Model	Colour	Size	Capacity
ESF9-0720S	Clear	Adult	1600 ml
ESF9-0721S	Clear	Child	500 ml
ESF9-0722S	Clear	Infant	280 ml

Clear silicon face masks

Masks can be used with every reanimation bag. Autoclavable

Model	Size
ESF9-0701	0
ESF9-0702	1
ESF9-0703	2
ESF9-0704	3
ESF9-0705	4
ESF9-0706	5



ESF9-0701



ESF9-0702



ESF9-0703



ESF9-0704



ESF9-0705



ESF9-0706

Silicon reanimation bags

100% latex free. Completely autoclavable. Made of pure, durable, silicon.
 All parts can be easily removed to ease cleaning or disinfecting.
 Complete with pressure release valve to minimize the risk of over inflating the lungs.
 The release valve can be locked if needed.



ESF9-0730S



ESF9-0731S



ESF9-0732S

Model	Colour	Size	Capacity
ESF9-0730S	Black	Adult	1600 ml
ESF9-0731S	Black	Child	500 ml
ESF9-0732S	Black	Infant	280 ml

Black silicon face masks

Masks can be used with every reanimation bag. Autoclavable

Model	Size
ESF9-0707	0
ESF9-0708	1
ESF9-0709	2
ESF9-0710	3
ESF9-0711	4
ESF9-0712	5



ESF9-0707



ESF9-0708



ESF9-0709



ESF9-0710



ESF9-0711



ESF9-0712

Silicon reanimation bag sets



ESF9-0720C5

ADULT Size

With autoclavable silicon face mask size 5 and reservoir, latex free



ESF9-0721C3

CHILD Size

With autoclavable silicon face mask size 3 and reservoir, latex free



ESF9-0722C1

INFANT Size

With autoclavable silicon face mask size 1 and reservoir, latex free



Silicon reanimation bag sets



ESF9-0720

ADULT Size

Reanimation bag 1600 ml complete with oxygen reservoir bag 2600 ml. Autoclavable, complete with 2 face masks size 3 and 5, 3 guedel airways size 1/3/4, oxygen tubing. Latex free



ESF9-0721

CHILD Size

Reanimation bag 500 ml complete with oxygen reservoir bag 2600 ml. Autoclavable, complete with 2 face masks size 2 and 3, 3 guedel airways size 0/2/4, oxygen tubing. Latex free



ESF9-0722

INFANT Size

Reanimation bag 280 ml complete with oxygen reservoir bag 600 ml. Autoclavable, complete with 2 face masks size 0 and 1, 3 guedel airways size 00/1/3, oxygen tubing. Latex free

Silicon reanimation bag sets with carrying case

Complete with: Autoclavable reanimation bag, oxygen reservoir bag with oxygen tubing, face mask, carrying case in nylon. Latex free



Model	Size
ESF9-0726	Set for adult
ESF9-0727	Set for child
ESF9-0728	Set for infant

Silicon reanimation bag set with carrying case

Complete with: Autoclavable silicon reanimation bag adult, silicon face masks sizes 0/1/2/3/4/5 and reservoir. Latex free.



ESF9-1050

Disposable PVC reanimation bag sets



ESF9-0734C5

ADULT Size

With single use mask size 5 and reservoir.



ESF9-0735C3

CHILD Size

With single use face mask size 3 and reservoir.



ESF9-0736C1

INFANT Size

With single use face mask size 1 and reservoir.

Disposable PVC resuscitating set with carrying case

Complete with:
PVC reanimation bag adult, single use masks sizes 3/5, guedel canula sizes 00/0/1/2/3/4 and reservoir, with comfortable bag for transportation. Latex free.



ESF9-1071

Double chamber reanimation bag

A reanimation bag with two chambers, completely made of autoclavable materials.

Model	Size
ESF9-2000R	Adult w/o mask 1500 ml
ESF9-2000M	Adult with mask 1500 ml
ESF9-2100R	Child w/o mask 500 ml
ESF9-2100M	Child with mask 500 ml





Silicon cushion masks

Silicone cushion masks are designed to use with manual resuscitators, automatic resuscitators and ventilators. Size 2, 3, 4, and 5 with a 22mm ISO connector. They can be autoclavable at 134 degree centigrade repeatedly. Transparent dome to observe patient's condition. Adjustable soft cuff inflation.

Mask dome: polysulphone.

Mask cuff: silicon rubber.

Mask connector: silicone rubber.

Hook: stainless steel.



ESF9-0772



ESF9-0773



ESF9-0774



ESF9-0775

Model	Size
ESF9-0772	22mm ISO connector 2
ESF9-0773	22mm ISO connector 3
ESF9-0774	22mm ISO connector 4
ESF9-0775	22mm ISO connector 5

CPR face mask

CPR face mask, patented one-way valve and directional diaphragm, allows no backflow of air or contaminants from patient to rescuer. Clear plastic, designed to fit adult or child. Soft inflated cuff. Mask not affected by extreme temperatures. Standard 22 mm. fittings - can be used with every kind of reanimation bag. Flexible dome folds into pocket size. Compact carrying case, with instructions.



ESF9-0719

Rebreathing bag

Rebreathing bag made of latex with loop, connection 22 mm, antistatic, autoclavable 121°C

References	Capacity
ESF9-0750	0,5 l
ESF9-0751	1 l
ESF9-0752	1,5 l
ESF9-0753	2 l
ESF9-0754	3 l
ESF9-0755	4 l
ESF9-0756	5 l
ESF9-0757	6 l



Anaesthesia face masks

Rubber face mask, black, inflatable pad with s.s. ring, antistatic, autoclavable 121°

Model	Size
ESF9-0760	0
ESF9-0761	1
ESF9-0762	2
ESF9-0763	3
ESF9-0764	4
ESF9-0765	5
ESF9-0766	6



ESF9-0760



ESF9-0761



ESF9-0762



ESF9-0763



ESF9-0764



ESF9-0765



ESF9-0766

Reservoir



072/RI4
Infant



072/RA4
Adult/Child

Empty carrying cases



ESF9-0727
With zip, handle and
transparent cover
32x12x11h cm



ESF9-0729
Transparent plastic case
32,5x12x10,5h cm



Medical bag with oxygen cylinder

Complete with: 0.5 L rechargeable steel oxygen bottle, fixed calibration reduction gear, autoclavable silicon reanimation bag, autoclavable silicon face mask, 3 Guedel canula, tongue puller, helicoidal mouth opener.
Size: 50x23x23 H cm

Model	Size
ESF9-1012	Set for adult
ESF9-1013	Set for child
ESF9-1014	Set for infant



Medical bag

Complete with: autoclavable silicon reanimation bag, autoclavable silicon face mask, 3 Guedel canula, tongue puller, helicoidal mouth opener.
Size: 40x16x21 H cm

Model	Size
ESF9-1018	Set for adult
ESF9-1019	Set for child
ESF9-1020	Set for infant



Intubation set for adult

Complete with:
1 Mc-intosh laryngoscope handle and blades sizes 1/2/3, 1 Magill forceps, 1Pean forceps, 1 scissor, 5 endo-tracheal tubes, 1 reamer for endotracheal tubes, 3 Guedel canula size 0/2/4, 1 mouth opener, 1 plaster roll, 5 suction catheters CH 14



ESF9-103050

Intubation set for child

Complete with:
Mc-Intosh laryngoscope handle child and blades sizes 0/1, operating scissor straight blunt cm. 10,5, hemostatic forceps straight Rochester-Pean cm 14, Magill forceps for child, mouth opener, Guedel canula size 0/1/2, endo-tracheal tubes uncuffed size 2.5/3/3.5/4, reamer, roll of wound plaster, 5 suction catheters CH 14



ESF9-103051

Resuscitating case

Resuscitating set in ABS container with fixing support, complete with: 0.5 L rechargeable steel oxygen bottle, fixed calibration reduction gear, autoclavable silicon reanimation bag, 3 autoclavable silicon face masks, 3 Guedel canula, tongue puller, helicoidal mouth opener. Dimensions 445x320x150 mm

Model	Size
ESF9-1006	Set for adult
ESF9-1007	Set for child
ESF9-1008	Set for infant



Reanimation kit

ABS case with fixing support, complete with: autoclavable silicon reanimation bag, oxygen reservoir with oxygen tubing, 2 silicon face masks, 3 Guedel canula sterile, mouth opener, tongue puller. Size 40x26x13 H cm

Model	Size
ESF9- REA0720	Set for adult
ESF9- REA0721	Set for child
ESF9- REA0722	Set for infant



Non rebreathing valves



072/PVW
60 cmH2O
For reanimation bag
adult, autoclavable



072/PVWC
40 cmH2O
For reanimation bag
child, autoclavable



072/PVW.01
60 cmH2O
For reanimation bag
adult, disposable



072/PVWC.01
40 cmH2O
For reanimation bag
child, disposable



Peep Valves and Diverters



072.P.10
PEEP VALVE
2-10 cmH₂O, ORANGE , 30M



072.P.20
PEEP VALVE
5-20 cmH₂O, BLUE, 30M



072.P.30
PEEP VALVE DIVERTER
Autoclavable. 30M/22F

Scissors

Multitask emergency scissors.
Made of s.s., plastic handle.
Sterilizable in autoclave.



Model	Size
ESF11-1257	15cm
ESF11-1258	18cm

Magill forceps

Made of Stainless Steel



Model	Size
ESF8-802	adult 25 cm
ESF8-801	child 20 cm
ESF8-800	infant 17 cm

Guedel canula

Canula to be used for resuscitation activities. Made of plastic material, sterile, disposable.

Model	Size	Model	Size
072/G000	000	072/G3	3
072/G00	00	072/G4	4
072/G0	0	072/G5	5
072/G1	1	072/G	Set of 6 pcs
072/G2	2		00/0/1/2/3/4



Isothermic blanket

Polyester isothermic emergency blanket.
Silver side protect patient against cold.
Gold side protect patient against heat.
Size: 150 x 220 cm



COP216

Mouth opener & tongue puller



PRA067
Mouth opener



PIN201
Tongue puller

Diagnostic penlight

Diagnostic lamp with pushbutton

Model	Size
ESF8-80003	without battery
ESF8-80004	with battery AAA 1,5 V type



Splinter removal pencil

Small metallic splinters eye tool remover.
One side with a magnetic point to remove any metallic splinter avoiding any risk for the eye.
On the other side an antiabrasive nylon loop for non metallic splinters



MAT022

Hot and Ice bags



ESF15-9050
Hot water bag 2000 ml



ESF15-9051
Ice bag cm25



Torniquet

Torniquet made of pure latex.
Lenght cm 33, box of 200 pcs



LAC001

Torniquet

Torniquet made in cotton with plastic hook.



CIN002

Sterile scalpel blades

Disposable carbon steel sterilized by ionizing radiation
Individually wrapped in aluminium foil and in boxes of 100 pcs



BL10



BL11



BL12



BL15



BL20



BL21



BL22



BL23



BL24

Scalpel handles for blades

Stainless steel.



ESF8-161

Size N°3 for blades 10/15



ESF8-162

Size N°4 for blades 20/24

Endo-tracheal tubes



Model	Description
INT313990	Uncuffed size 2.5
INT313992	Uncuffed size 3.0
INT313994	Uncuffed size 3.5
INT313996	Uncuffed size 4.0
INT313998	Uncuffed size 4.5

Model	Description
INT314000	Cuffed size 5.0
INT314002	Cuffed size 5.5
INT314004	Cuffed size 6.0
INT314006	Cuffed size 6.5
INT314008	Cuffed size 7.0
INT314010	Cuffed size 7.5
INT314012	Cuffed size 8.0
INT314014	Cuffed size 8.5

Intubation reamer




INT314024

Venturi Mask

With 6 Venturi diluters, 2.10 m O₂-sure flow tubing and 15 cm flex tubing

Model	Size
ESF5-1710	Adult
ESF5-1711	Paediatric





Oxygen mask

Oxygen mask with 2.10 m O₂-sure flow tubing


Model	Size
ESF5-1041	Adult
ESF5-1043	Paediatric
ESF5-1040	Adult without tube
ESF5-1042	Paediatric without tube



Nasal cannula

With 2.10 m O₂-sure flow tubing, adult

Model	Size
ESF5-1012	Adult
ESF5-1013	Paediatric



High concentration mask

Oxygen mask, non-rebreathing, with bag, check valve, safety vents closed and 2.10 m O₂-sure flow tubing

Model	Size
ESF5-1095	Adult
ESF5-1096	Paediatric



Aerosol mask

Complete set of aerosol mask, nebulizer bulb and O₂-sure flow tubing

Model	Size
ESF5-5010	Adult set
ESF5-5009	Paediatric set
ESF5-5011	Nebulizer Bulb
ESF5-5012	Adult mask
ESF5-5013	Paediatric mask
ESF5-5015	Oxygen sure flow tubing 2.10 m



Aerosol mask adult size for tracheostomized patients

Tracheostomy mask for 22mm corrugated tubing, with elastic strap and 360° connector



ESF5-5090

Portable cylinders for oxygen therapy

The portable cylinders for oxygen therapy are F.U. compliant. Sold empty. Made in steel, with valve E2000 and connection E17. Refillable and including a test certificate.
Optional: protection case

Model	Size
ESF9-15029	0,5 L
ESF9-15031	1 L
ESF9-15032	2 L
ESF9-15033	3 L
ESF9-15035	5 L
ESF9-15037	7 L
ESF9-15038	10 L



Reanimation set

Complete with:
rechargeable oxygen steel bottle 5L , reducer with flowmeter and scrubber, silicon black reanimation bag , face mask size 3 and oxygen tube. Mounted on a steel trolley with castors

Model	Size and Weight
ESF9-3003	35x17x68H cm 11,5 kg





Piston pressure reducer with flowmeter 0-15 l/min and humidifier

Piston pressure reducer with fixed calibration, equipped with humidifier and flowmeter.
 Delivery flow max 4bar.
 Pressure gauge Ø 42, scale 0-315 bar;
 Precision class 2.5.
 Cylinder connection UNI4406.
 Output connection with 6mm diameter hose connector.



ESF9-12214

Piston pressure reducer with flowmeter 0-15 l/min

Piston pressure reducer with fixed calibration, equipped with flowmeter.
 Delivery flow max 4bar.
 Pressure gauge Ø 42, scale 0-315 bar;
 Precision class 2.5.
 Cylinder connection UNI4406.
 Output connection with 6mm diameter hose connector



ESF9-12215

O2 regulator

Adjustable output flow control 0-12 L/min.
 Pressure gauge Ø42,
 Scale 0-315bar;
 Precision class 2.5
 Cylinder connection UNI4406



ESF9-12201

Pressure reducer

Piston pressure reducer with fixed calibration.
 To be used with cylinders equipped with medical gas reduction handwheel.
 Chrome brass main body with safety valve.
 Delivery flow: 6 L/min.
 Pressure gauge:
 diameter 42mm, scale 315bar, precision class 2.5.
 Output connection with 6mm diameter hose connector



ESF9-12246



ESSE3 srl, Via Garibaldi 30
 14022 Castelnuovo D.B. (AT)
 Tel +39 011 99 27 706
 Fax +39 011 99 27 506
 e-mail esse3@chierinet.it
 web: www.esse3-medical.com

