



DATA SHEET

Models: 21010 ... 21023

SHYFY 'gurgical gowns'

Specialistic, water repellent, breathable gowns for a barrier effect to fluids and bacteria

- cotton cuffs 8 cm
- ultrasonic welding
- neck velcro closure
- 4 waist straps for internal and external closure

Pre folded


For easy programmed dressing
Welded reinforcement (only 21015, 21016)

Water-proof, breathable reinforcement on front and on sleeves ultrasonic welded.


All gowns include two dry-paper towels.
Conforms to UNI EN ISO 13795.



Cod.	Description	
	SURGICAL GOWNS 35 g	
21010	Size M - length 115 cm	box 30
21011	Size L - length 130 cm	box 30
21012	Size XL - length 135 cm	box 30
	SURGICAL GOWNS 50 g	
21020	Size M - length 120 cm	box 30
21021	Size L - length 130 cm	box 30
21022	Size XL - length 140 cm	box 30
21023	Size XXL - length 150 cm	box 30
	REINFORCED GOWNS 35 g	
21015	Size L - length 130 cm	box 30
21016	Size XL - length 135 cm	box 30



SURGICAL GOWNS 35/50 g



REINFORCED GOWNS 35 g



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

