

ESF LED THIN X-RAY FILM VIEWERS NEGATOSCOPE LED ULTRA MINCE



- Frame made of aluminum alloy and plastic. Viewing area 36x44 cm
- Wall-mounting
- Cadre en alliage d'aluminium et plastique. Zone de visualisation 36x44 cm
- Montage mural

CODE	DIMENSIONS
ESF2-041/00	47 x 54,5 x 2,4 cm
ESF2-081/00	83 x 54,5 x 2,4 cm



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail: danielec@esse3-medical.com
web: www.esse3-medical.com

