

ESF WARD ROOM ACCESSORIES



ESF2-694 I.V. STAND ON CASTORS

Adjustable height from 115 to 220 cm.

STAINLESS STEEL BASE: Foldable base with 5 castors

- /60 S.S. tube, 2 plastic hooks,
- /61 S.S. tube, 4 S.S. hooks
- /75 S.S. tube, 4 S.S. hooks, antistatic
- /80 S.S. tube, 4 S.S. hooks.
- /81 S.S. tube, 4 S.S. hooks, antistatic castors

PLASTIC BASE : Base with 5 castors

- /11 S.S. tube, 4 S.S. hooks
- /12 S.S. tube, 2 plastic hooks
- /13 S.S. tube, 4 plastic hooks
- /13A S.S. tube, 4 S.S. hooks

ESF2-694 PIED A SERUM MOBILE

Hauteur réglable de 115 à 220 cm.

PIÈTAMENT INOXYDABLE Piètement pliant avec 5 roulettes

- /60 Tubeinox, 2 crochets en plastique
- /61 Tube inox, 4 crochets en acier inox,
- /75 Tube inox, 2 crochets en acier inox, roulette A/S
- /80 Tube inox, 4 crochets en acier inox.
- /85 Tube inox, 4 crochets en acier inox, roulette A/S

PIÈTAMENT PLASTIQUE Piètement avec 5 roulettes

- /11 Tube inox, 4 crochets en acier inox
- /12 Tube inox, 2 crochets en plastique
- /13 Tube inox 4 crochets en plastique
- /13A Tube inox, 4 crochet en acier inox



SOILED LINEN TROLLEY

Chrome plated frame.

Mounted on Ø 80 mm casters, removable bag.

Knock down construction

CHARIOT POUR LINGE SALES

Structure en acier chromé.

Monté sur roulettes Ø 80 mm, sac amovible.

Structure démontable

CODE	
ESF2-840/00	Single storage cart
ESF2-841/00	Double storage cart



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail: danielec@esse3-medical.com
web: www.esse3-medical.com



ESF WARD ROOM ACCESSORIES



MOUNTING STEP

Chrome plated steel frame with plastic feet.
Available also in stainless steel .
Step covered with anti slip rubber. Knock down construction.

CODE	
ESF2-850/01	S/S. 1 step
ESF2-851/01	S/S. 2 steps

MARCHEPIED

Structure en tube d'acier chromé.
Disponible aussi en version acier inoxydable
Marche avec revêtement antidérapant. Structure démontable



WARD TROLLEY / CHARIOT DE SERVICE

CODE	SIZE/DIMENSIONS
ESF2-1050/01	1 drawer, ss bowl, 70x50x80H cm
ESF2-1051/01	2 drawers, ss bowl, 70x50x80H cm

INSTRUMENT TROLLEY

- Made of stainless steel 18/10.
- Complete with: drawer and S.S. bowl 18/10 ø 30 cm.
- Dimensions: 50 x 70 x 80cm

CHARIOT A INSTRUMENT

- Structure en acier inox 18/10.
- Complet avec: tiroir et cuvette inox ø 30cm.
- Dimensions: 50 x 70 x 80cm