

# ESF GYNECOLOGICAL COUCH



## **ESF2-360/10 GYNAECOLOGICAL COUCH**

Constructed from Chrome-plated round steel tube.

Adjustable back-rest and leg sections.

Top is upholstered and covered with washable plastic material flame retardant.

Complete with S.S. steel bowl  $\varnothing$  32 cm., and two standard padded leg holders.

DIMENSIONS: 180x55x80 H cm.

Knock down construction.

## **ESF2-360/10 DIVAN GYNECOLOGIQUE**

Structure en acier chromé.

Têteière et plateau de jambes réglables.

Matelas rembourré avec revêtement anti-feu lavable.

Livré avec: 1 cuvette en acier inox  $\varnothing$  32 cm. et appui-jambes standard rembourré.

DIMENSIONS: 180x55x80 H.cm

Structure démontable.

# ESF GYNECOLOGICAL COUCH



## GYNAECOLOGICAL COUCH

Constructed from chrome-plated rectangular steel tube.

Adjustable back-rest and leg sections.

Top is upholstered and covered with washable plastic material flame retardant.

Complete with S.S. steel bowl  $\varnothing$  32 cm.

DIMENSIONS: 180x55x80 H cm.

Knock down construction.

**ESF2-359/00** knee crutches

**ESF2-359/01** standard padded leg holders

## DIVAN GYNECOLOGIQUE

Construction en acier chromé de section carrée.

Tête et plateau de jambes réglables.

Matelas rembourré avec revêtement anti-feu lavable.

Livré avec: 1 cuvette en acier inox  $\varnothing$  32 cm.

DIMENSIONS: 180x55x80 H.cm

Structure démontable.

**ESF2-359/00** appui-jambes anatomiques

**ESF2-359/01** appui-jambes standard



ESSE 3 Via Garibaldi 30  
14022Castelnuovo D.B. (AT)  
tel +39 011 99 27 706  
fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web : [www.4ci.com](http://www.4ci.com)

