



ESF2-2863/00

Delivery bed



MEDICAL FURNITURE

ESF2-2863/00

Delivery bed



Cod. **ESF2-2863/00**



Optional:



The birthing bed has been designed focusing particularly on ergonomics, safety and mobility.

Suitable for all stages of labour and for hospitalization pre and after labor. Multiple delivery positions in a single unit result in a great level of comfort for the mother and medical personnel. The midwife can easily access the mother from different positions thus improving ergonomics and safety. Adjustable leg rests for ergonomic positioning with locking of the clamp into different positions. Bed ends removable instantly via quick release mechanism. Wide mattress base. More comfortable anatomical leg rest. Side handles for enhanced pushing force. Side rails foldable. Easy to get in and out thanks to low minimum height.

It can also be used for transportation within the hospital. Easy-to-clean structure. Hydraulic birthing bed is designed to ensure the best possible usage for the product without the need for electricity power supply.

All adjustments in this model are manual.

The height of the bed can be adjusted with a hydraulic pedal, the back section adjustment is with crank infinitely adjustable located at bed head for assisted elevation and the leg section slides manually in and out. Trendelenburg on the seat has been realized using 2 racks to assure safety and stability during labour operation.

These features provides enhanced possibilities to use various positions during labour. Electric birthing bed offer a very silent motor for electrical variable height and different patient positions are adjustable with gas pumps and the leg section slides manually in and out.

Trendelenburg on the seat has been realized using 2 racks to assure safety and stability during labour operation. These features provides enhanced possibilities to use various positions during labour.

Technical Features:

- Width 1000 mm
- Mattress base (width) 800 mm
- Height adjustment 600 - 950 mm
- Length 1700 - 2200 mm
- Safe working load (SWL) 230 kg
- Trendelenburg 12 °
- Back section adjustment 0 - 80 °
- Castors size 150 mm
- Mattress base (sections) 3 pcs



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

