

ESF DONNOR BEDS



DONNOR / TREATMENT COUCH WITH VARIABLE HEIGHT

Adjustable head and leg sections in any position.
Top is up holstered and covered with washable plastic material in white color.
Frame made in sheet steel epoxy painted with adjustable feet.

Variable height from 56 to 93 cm.

Armrests re movable

ESF2-538/00 variable height by electric motor

ESF2-538/10 variable height by hydraulic pump



DIVAN MULTIFONCTIONNEL HAUTEUR VARIABLE

Rèleve buste et plateau des jambes réglables.

Matelas rembourré avec revêtement lavable en couleur blanc.

Construction en acier émaillé.

Hauteur réglable de 56 à 93 cm.

Appuie-bras amovibles.

ESF2-538/00 hauteur réglable avec moteur électrique

ESF2-538/10 hauteur variable par pompe hydraulique

ESF DONNOR BEDS



BLOOD DONOR CHAIR

Chrome plated steel tube.

Leg, seat, back and arm rest padded and covered with washable plastic material.

Simultaneous movement of back and leg rest.

Adjustable arm rests.

Options: 1) backrest height adjustable; 2) backrest height and leg rest adjustable; 3) backrest height, leg rest and arm rest adjustable.

FAUTEUIL REGLABLE POURPRISE DE SANG

Structure en acier chromé, entièrement revêtu de plastique lavable.

Mouvement simultané du dos et des jambes, appuie-bras réglable

) Fermée 100x55x120H cm / ouverte 175x55x50H cm



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.4ci.com

