



# MULTIFUNCTIONAL CART



**Cod. 816191/01**

Multipurpose cart with : 3 drawers in painted steel, 1 Drawer in painted steel, 1 ABS extending table, 1 integrated tray with raised rim on left side.

The self-bearing monocoque body is made in steel sheet coated with anti-mould epoxy powders.

The left compartment features four trays with metal rim.

Drawers are made in steel, with anti-mould epoxy powder coating, and are equipped with: Central locking system, Silent sliding system on metal tracks and rollers, End stoppers.

The tabletop is made in technopolymers, and has no sharp edges for ease of use cleaning and disinfecting.

Push handle in stainless steel, fully integrated in the frame.

4 casters, of which 2 with brake, Ø125mm (anti-static on demand), non-marking, anti-thread, mounted on swiveling bearings, 4 PVC non-marking bumpers.

Optional:

**DR005** - Holding frame with film viewer  
40 x 43 cm

## *Stainless Steel*



80 X 62 X 100 H (cm)



**ESSE3** srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)





# CHARIOT MULTIFONCTIONNEL



**Cod. 816191/01**

Chariot multi-usage avec: 3 tiroirs en acier laqué, 1 tiroir en acier laqué, 1 plateau extractible en ABS, 1 plan porte-objets avec protection avant sur le côté gauche. Le cadre autoporteur est réalisé en tôle d'acier laqué à l'aide de poudres époxy anti-moisissure.

Dans le compartiment gauche sont logés quatre petits plans avec une protection avant métallique.

Les tiroirs en acier, laqués à l'aide de poudres époxy anti-moisissure, sont livrés avec : Système de fermeture avec clé centrale, Système de coulissement silencieux sur glissières métalliques et billes, Butée de fin de course.

Le plan de travail en technopolymères ne présente aucune arête pour faciliter le nettoyage et la désinfection.

Poignée réalisée en acier inox parfaitement intégrée dans la structure.

Avec base 4 roues, dont 2 avec frein, de 125 mm II (antistatiques sur demande), anti-traces, ne s'effilochant pas, montées sur roulements pivotants, 4 pare-chocs en PVC anti-rayure.

En option :

**DR005** - Rehausse avec négatoscope  
40 x 43 cm

## *Acier Inoxidable*



80 X 62 X 100 H (cm)



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

