



# MULTIFUNCTIONAL CART



## Cod. **ESF1-6192**

Multifunction cart with : 5 Drawers in painted steel, 1 ABS extending table, 1 white, ABS, 10 liter waste bin, 1 red ABS catheter and probe Holder The self-bearing monocoque body is made in steel sheet coated with anti-mould epoxy powders. The left compartment features four trays with metal rim. Drawers are made in steel, with anti-mould epoxy powder coating, and are equipped with : Central locking system, Silent sliding system on metal tracks and rollers, End stoppers. 1 drawer mm. 615 x 470 x 68 H, 4 drawers mm. 615 x 470 x 147 H . The tabletop is made in technopolymers, and has no sharp edges for ease of use when cleaning and disinfecting. Push handle in stainless steel, fully integrated in the frame 4 casters, of which 2 with brake, Ø 125mm (anti-static on demand), non-marking, anti-thread, mounted on swiveling bearings, 4 PVC non-marking bumpers.



80 X 62 X 100 H (cm)



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)





# CHARIOT MULTIFONCTIONNEL



## Cod. **ESF1-6192**

Chariot multi-usage avec : 5 tiroirs en acier laqué, 1 plateau extractible en ABS, 1 poubelle en ABS de 10 litres blanche, 1 porte-cathéter et sondes en ABS rouge. Le cadre autoporteur est réalisé en tôle d'acier laqué à l'aide de poudres époxy anti-moisissure. Dans le compartiment gauche sont logés quatre petits plans avec une protection avant métallique. Les tiroirs en acier, laqués à l'aide de poudres époxy anti-moisissure, sont livrés avec : Système de fermeture avec clé centrale, Système de coulissement silencieux sur glissières métalliques et billes, Butée de fin de course 1 tiroir de 615 x 470 x 68 mm/ H, 4 tiroirs de 615 x 470 x 147 mm/ H. Le plan de travail en technopolymères ne présente aucune arête pour faciliter le nettoyage et la désinfection. Poignée réalisée en acier inox parfaitement intégrée dans la structure. Avec base 4 roues, dont 2 avec frein, de 125 mm Ø (antistatiques sur demande), anti-traces, ne s'effilochant pas, montées sur roulements pivotants, 4 pare-chocs en PVC anti-rayure.



80 X 62 X 100 H (cm)



**ESSE3** srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

