



# MEDICINE TROLLEY



## Cod. **ESKR-MT/720**

- ABS body
- Metal ball-bearing rail.
- Dust basket
- Utility container
- File bag
- Centralized lock
- Inner with partitions can be well organized
- Luxurious Noiseless casters with brake
- Powder coated column



630x475x920h (mm)



**ESSE3** srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

