



EG6!\$&!\$' `VUF]UHf]W

9`YWF]W6 YX'

k jh ` : jj Y' : i bW]cbg



MEDICAL FURNITURE



ESB-02-03 bariatric

Electric Bed with Five Functions



Technical specifications:

- Full Length : 223 cm
- Full Width : (including the side rails 105 cm) (Without the side rails 99 cm)
- Bed base height range : (Maximum height 70 cm) (Minimum height 45 cm)
- Patient touching surface : 90 x 196 cm
- Bed base angle position Adjustable :
- Back section panel (0 ~ 75°)
- Leg section panel (0 ~ 75°)
- Trendelenburg (0 ~ 10°)
- Anti-Trendelenburg (0 ~ 10°)
- Maximum safe load : 250 kg
- Diameter of Centrally Silent Castor : 12,5 cm
- Patient controls in side rails
- Instant CPR release lever
- Standard telescopic IV pole with 4 hooks