



ESF1-300/VH

4 Section bed with variable height



MEDICAL FURNITURE



ESF1-300/VH

Technical Specifications:

FOUR SECTIONS FOWLER BED, WITH ALL WELDED ENAMELLED MAIN FRAME MADE OF RECTANGULAR SECTION STEEL TUBE.

REMOVABLE CHROME PLATED BED ENDS FITTED WITH 10 MM. LAMINATED PLASTIC PANELS.

PLATFORM FITTED WITH EPOXY-PAINTED STEEL WELDMESH IN FOUR SECTIONS.

ADJUSTABLE BACK-REST AND LEG-SECTION BY MEANS OF TWO INDEPENDENT CRANCKS WITH RETRACTABLE HANDLE.

VARIABLE HEIGHT ADJUSTABLE BY HYDRAULIC PUMP.

BUMPERS AT ALL FOUR CORNERS.

MOUNTED ON 4 Ø125MM CASTORS, TWO WITH BRAKES.

KNOCK DOWN CONSTRUCTION. 50/75CM

Length: 215.00CM

Width: 92.00CM

Height: 50.00CM

Accessory:

MATTRESS

