

ESF OPERATING TROLLEY

Constructed from stainless steel 18/10.
Shelves with up stands on four sides.
Push handle.
Mounted on \varnothing 80 mm. casters.
Knock down construction

Structure en acier inox 18/10.
Deux plateaux avec 4 bords relevés.
Poi gnée.
Monté sur roulettes \varnothing 80 mm.
Structure démontable.



MEDICAL TROLLEY / TABLE DE PANSEMENT

| CODE | SIZE/DIMENSIONS |
|--------------------|-----------------|
| ESF2-869/00 | 70x50x80H cm |



INTRUMENT TROLLEY / GUERIDON PORTE INSTRUMENTS

| CODE | SIZE/DIMENSIONS |
|--------------------|-----------------|
| ESF2-866/00 | 60x40x80H cm |
| ESF2-867/00 | 70x50x80H cm |
| ESF2-868/00 | 90x60x80H cm |



DRESSING TROLLEY / TABLE DE PANSEMENT

| CODE | SIZE/DIMENSIONS |
|---------------------|----------------------------------|
| ESF2-1050/00 | 1 drawer, ss bowl, 70x50x80H cm |
| ESF2-1051/00 | 2 drawers, ss bowl, 70x50x80H cm |



DRESSING TROLLEY / TABLE DE PANSEMENT

| CODE | SIZE/DIMENSIONS |
|---------------------|---------------------------------------|
| ESF2-1050/00 | Rack, wastebin, ss bowl, 90x60x80H cm |

ESF OPERATING ACCESSORIES



MAYO TABLE / TABLE DE MAYO

| CODE | SIZE/DIMENSIONS |
|--------------------|--|
| ESF2-481/00 | Height adjustment 74-104 cm Tray size 64 x 43 x 3 cm - Hauteur réglable de 80 à 115 cm. Dimensions du plateau 64 x 43 x 3 cm |



KICK BUCKET / BAC A DECHETS

| CODE | SIZE/DIMENSIONS |
|--------------------|---|
| ESF2-541/00 | Manufactured from 18/10 stainless steel. Mounted on 5 casters Ø 40 mm. - Acier inoxydable 18/10. Monté sur 5 roulettes Ø 40 mm. |



PLASTER BOWL / CUPULE POUR PLATRE

| CODE | SIZE/DIMENSIONS |
|--------------------|---|
| ESF2-551/00 | Manufactured from 18/10 stainless steel w/o welding.. - Acier inoxydable 18/10 sans soudure. |



BOWL / CUVETTE

| CODE | SIZE/DIMENSIONS |
|--------------------|---|
| ESF2-535/00 | Stainless steel 18/10 Ø 32 cm. Capacity 4 lt. - Acier inoxydable 18/10, Ø 32 cm. Capacité 4 lt. |

ESF OPERATING ACCESSORIES



BOWL'S STAND

Constructed from chrome plated steel tube.
Variable height.
Stainless Steel bowl.

PORTE CUPULES

Colonne de support en tube d'acier chromé.
Hauteur réglable.
Cuvette en acier inoxydable

| CODE | | |
|-------------|---------------------------------------|--|
| ESF2-522/00 | Single basin stainless , steel base | Simple bassin , base en acier inoxydable |
| ESF2-523/00 | Single basin , plastic base | Simple bassin , base en plastique |
| ESF2-532/00 | double basin stainless , steel base | Double bassin , base en acier inoxydable |
| ESF2-533/00 | double basin stainless , plastic base | Double bassin , base en plastique |



| CODE | WASTE BIN Constructed from stainless steel 18/10. Removable inner container of plastic material. | POUBELLE Acier inoxydable 18/10. Seau amovible en plastique. |
|-------------|---|--|
| ESF2-524/00 | Capacity 20 lt | Capacité de 20 l |
| ESF2-545/00 | Capacity 14 lt | Capacité de 14 l |
| ESF2-526/00 | Capacity 5 lt | Capacité de 5 l |

ESF OPERATING ACCESSORIES



ESF2-490/00 ANAESTHESIA TROLLEY

Constructed from stainless steel 18/10. Upper section with a large open space. Six drawers and storage compartment complete with shelf and locking door. Mounted on \varnothing 80 mm. casters. DIMENSIONS : 40x60x110 (H) cm.

ESF2-490/00 GUERIDON D'ANESTHESIE

Structure en acier inoxydable 18/10 comprenant : 1 étagère, 6 tiroirs, 1 armoire de rangement avec serrure. Monté sur roulettes diam. 80 mm. DIMENSIONS: 40x60x110 (H) cm.



WARD SCREEN

Chrome-plated steel tubular frame, plastic curtains washable and flame retardant. Mounted on casters. Knock down construction.

PARAVENT

Tube d'acier chromé panneaux anti-feu et lavables. Mobile sur roulettes. Construction démontable

| CODE | |
|--------------------|--------------------------------------|
| ESF2-824/00 | 2 panels , size 100x42x167 cm |
| ESF2-825/00 | 3 panels , size 150x42x167 cm |
| ESF2-826/00 | 4 panels , size 200x42x167 cm |



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.4ci.com

