



HYDRAULIC MAYO TABLES



MEDICAL FURNITURE



HYDRAULIC MAYO TABLES



Cod. **ESF2-482/00**

- Stainless steel construction
- Variable height with gas pump
- 4 castor base
- Max. load: 60 kg
- Swivel casters: 50 mm
- Dimensions: 45 x 67,5 x 90/125 H. cm

Cod. **ESF2-485/00**

- Stainless steel construction
- Variable height with gas pump
- 4 castor base
- Max. load: 60 kg
- Swivel casters: 50 mm
- Dimensions: 45 x 67,5 x 90/125 H. cm



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

