



# ESX-01 R

## Baby Cot




MEDICAL FURNITURE




# EGL!\$%&


.....6Uvm7ch



Trendelemburg -12°



Anti-Trendelemburg +12°



Height from vertical post to the edge 600 mm



Unlinable mattress



Removable Acrylic bassinets



3" PVC covered casters

## Technical Features

<b>Frame</b>	Epoxy powder coated steel frame, round steel, 1.5mm thickness. Tube diameter: 25mm
<b>Crib bassinet</b>	Transparent Perforated Acrylic bassinet with name card holder
<b>Loading weight</b>	30kg
<b>Tr./Antitr.</b>	Easily adjustable by screws, which allows tilt angle of +/- 12°
<b>Casters</b>	4 PVC covered casters, 2 with brake, diameter 76.2mm
<b>Mattress cover</b>	Waterproof and urine-proof mattress PVC cover, easily cleanable by disinfectants
<b>Mattress</b>	Mattress core is brand new high density polyurethane foam, 40kg/m3. The unlined mattress is washable and sterilizable. 625x375x45mm(LxWxH)
<b>Size</b>	900x525x650-1000mm



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

