



ES-OXY-20SA

Oxygen Concentrator



OPERATING THEATRE




Application

- Medical use : connect with ventilator and anaesthesia machine.
- Connecting small central oxygen supply.
- Oxygen supply to clinics and small hospitals with 5-50 beds.
- Optional Booster Pump: 2-8 Bar.

Technical Specification

Oxygen Flow: 0-20LPM
 Oxygen Concentration: 93% (±3%)
 Outlet Pressure: 0.04-0.07Mpa (6-10 PSI)
 Power: AC220V±10%,50/60HZ
 Sound Level: <60db
 Power Consumption: < 1100w
 Net Weight: 95Kgs
 Size: 680x500x1580

Schematic



ESSE3 srl, Via Garibaldi 30
 14022 Castelnuovo D.B. (AT)
 Tel +39 011 99 27 706
 Fax +39 011 99 27 506
 e-mail esse3@chierinet.it
 web: www.esse3.dreamgest.com

