



# HIGH PRESSURE OXYGEN GENERATORS

## Special for surgery



### Cod. ES-OXY10HF

Basic Functions: Switch, Oxygen gauge, Pressure regulating valve, LCD Display, Sound alarm, Power lamp, Temperature lamp, High/low pressure lamp, High/low oxygen lamp

#### High Pressure Oxygen Generator 0.14 Mpa fixed

- Applied in surgical especially connecting with ventilator or anesthesia machine
- Automatic Dehumidify
- Features: Automatic on and off
- Manually remove alarm
- Emergency stop button
- Audio and visual alarm for high & low pressure
- Full digital control system

#### Technical parameters:

- Power consumption: 650W
- Oxygen Flow(L/min): 10
- Air Flow(L/min): 10
- Oxygen Outlet Pressure(Mpa): 0.14Mpa
- Air outlet pressure(Mpa): 0.14Mpa
- Oxygen Concentration:  $\geq 90\%$
- Sound Level:  $\leq 60\text{db}$
- Net Weight: 40KG





# HIGH PRESSURE OXYGEN GENERATORS

## Special for surgery

### Cod. ES-OXY10HC

Basic Functions: Switch, Oxygen gauge, Pressure regulating valve, LCD Display, Sound alarm, Power lamp, Temperature lamp, High/low pressure lamp, High/low oxygen lamp

#### High Pressure Oxygen Generator 0.14 ~ 0.4 Mpa adjustable

- Applied in surgical especially connecting with ventilator or anesthesia machine
- Automatic Dehumidify
- Features : Automatic on and off
- Manually remove the alarm
- Emergency stop button
- Audio and visual alarm for high & low pressure
- Full digital control system

#### Technical parameters:

- Power consumption:850W
- Oxygen Flow(L/min):10
- Air Flow(L/min):10
- Oxygen Outlet Pressure(Mpa): 0.14-0.4Mpa
- Air outlet pressure(Mpa): 0.14-0.4Mpa
- Oxygen Concentration:  $\geq 90\%$
- Sound Level:  $\leq 60\text{db}$
- Net Weight: 65KG
- Optional: 15L & 20L



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

