



ES-SM/229A

Aspirator for Plaster Saw



OPERATING THEATRE



Versatile and compact vacuum cleaner equipped with a special filter for fine dust of plaster and synthetic material.
 The powerful 1600 W turbine generates an air flow of 61 liters / sec. and a depression of 244 mbar.
 The special filter blocks 99.97% of particles with a diameter greater than 0.5 µm.
 The automatic start and stop device is controlled directly by the switch located on the saw.
 Power supply: 220-240 Volt 50 / 60Hz - Extractor capacity: 15 liters
 Quantity of extractor air: 61 l / s - Exhaust noise level: 95 dB (A).

ACCESSORIES

- ES-SM/229M Plaster saw stand complete with height-adjustable handle and support (box 1 pc)
- ES-SM/229Z Paper bag for aspirator (box 5 pc)
- ES-SM/229O Fell bag for aspirator (box 5 pc)
- ES-SM/229W HEPA filter for aspirator (box 1 pc)



ESSE3 srl, Via Garibaldi 30
 14022 Castelnuovo D.B. (AT)
 Tel +39 011 99 27 706
 Fax +39 011 99 27 506
 e-mail esse3@chierinet.it
 web: www.esse3.dreamgest.com

