



# ESS-MAESTRO

## Bipolar TURP with Seal



OPERATING THEATRE

PROGRAM  
Program



## Quick Reference:

ESS-MAESTRO is a micro controller based model with Touch screen Display. Full power output at low cost. This is the right model for the customers who required 400w micro controller model.

Power Digital: 400 W

Programs: 100

Capability: All surgeries Top of the line model.

## Technical Specification:

Max Power: 400W.

Interface: Graphical LCD Display

Programs: 100

Monopolar Cut Modes: Pure Cut, Blend1, Blend2, Endo Cut.

Monopolar Coagulation Modes: Spray, Dessicate, Fulgurate.

Bipolar Modes: Bi-Cut ,Bi-Coag

TURP Modes: Bi-TURP-Cut, -TURP-Coag


Sealer mode: Seal

Tissue Response Circuit: Present.

Patient Plate Monitoring System: Present.

Dimension (WxHxD) & Weight: 300x150x340 mm / 5.5 kg

## Standard Accessories:



Monopolar Double Pedal Footswitch  
/ Bipolar Single Pedal Footswitch



Patient Plate With Cable  
Cord Autoclavable




Silicon Patient Plate (Optional)



Chuck Handle With Electrodes  
Autoclavable



Bipolar Forcep Autoclavable  
With Cable



Disposable Hand Switch Pencil



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3.dreamgest.com](http://www.esse3.dreamgest.com)

