



## SURGICAL LED OPERATING LIGHTS



### >> ES-HEX/CT5

Adjustable Colour Temperature

#### Special Features

Shadowless Light, Memory Function, High quality sterilizable handle, With protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.

#### Features

-  Memory Function
-  Endo Mode



#### Optional

-  Wireless Wall Controller
-  Battery Backup

#### Technical Specifications

Intensity	140,000 ±10%
Size of Light Field	12-30 cm
Colour Temperature	4000-5000 K
Colour Rendering Index	93 RA
LED Life	>50000 Hrs.
Number of LED	95
Diameter of Light	500mm
Brightness Control	Capacitive Touch Panel
Power Supply	220V/50Hz AC
Colour of LED	White & Yellow (Double Colour)
Focusing	Adjustable



**Mobile Model**

**Mobile Model**  
with battery backup



### >> ES-HEX/CT4

Adjustable Colour Temperature

#### Special Features

Shadowless Light, Memory Function, High quality sterilizable handle, With protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.

#### Features

-  Memory Function
-  Endo Mode



#### Optional

-  Wireless Wall Controller
-  Battery Backup

#### Technical Specifications

Intensity	110,000 ±10%
Size of Light Field	12-30 cm
Colour Temperature	4000-5000 K
Colour Rendering Index	93 RA
LED Life	>50000 Hrs.
Number of LED	76
Diameter of Light	500mm
Brightness Control	Capacitive Touch Panel
Power Supply	220V/50Hz AC
Colour of LED	White & Yellow (Double Colour)
Focusing	Adjustable



**Mobile Model**

**Mobile Model**  
with battery backup



# SURGICAL LED OPERATING LIGHTS

## >> ES-HEX/CT3

Adjustable Colour Temperature

### Special Features

Shadowless Light, Memory Function, High quality sterilizable handle, With protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.

### Features

- Memory Function
- Endo Mode



Mobile Model

Mobile Model with battery backup

### Optional

- Wireless Wall Controller
- Battery Backup

### Technical Specifications

Intensity	80.000 ±10%
Size of Light Field	12-30 cm
Colour Temperature	4000-5000 K
Colour Rendering Index	93 RA
LED Life	>50000 Hrs.
Number of LED	57
Diameter of Light	500mm
Brightness Control	Capacitive Touch Panel
Power Supply	220V/50Hz AC
Colour of LED	White & Yellow (Double Colour)
Focusing	Adjustable

