



LED SURGICAL LIGHTS

ES-NOVA6

Adjustable Colour Temperature

Special Features

Shadowless Light, Memory Function, High quality sterilizable handle. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.



M Memory Function

E Endo Mode

Technical Data

Intensity	160,000 Lux \pm 10%
Size of Light Field	12-30 cm
Colour Temperature	4000-5000 K
Colour Rendering Index	93 RA
LED Life	>50000 Hrs.
Number of LED	126
Diameter of Light	600mm
Brightness Control	Capacitive Touch Panel
Power Supply	220V/50Hz AC
Colour of LED	White & Yellow (Double Colour)
Focusing	Adjustable
Battery Backup	Optional



Mobile Model

Mobile Model with battery backup

ES-NOVA5

Adjustable Colour Temperature

Special Features

Shadowless Light, Memory Function, High quality sterilizable handle. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.



M Memory Function

E Endo Mode

Technical Data

Intensity	140,000 Lux \pm 10%
Size of Light Field	12-30 cm
Colour Temperature	4000-5000 K
Colour Rendering Index	93 RA
LED Life	>50000 Hrs.
Number of LED	105
Diameter of Light	600mm
Brightness Control	Capacitive Touch Panel
Power Supply	220V/50Hz AC
Colour of LED	White & Yellow (Double Colour)
Focusing	Adjustable
Battery Backup	Optional



Mobile Model

Mobile Model with battery backup



ESSE3 srl, Via Garibaldi 30
14022 Castelnovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail: danielec@esse3-medical.com
web: www.esse3-medical.com

