



# ESL-50 & ESL-60

## Emergency Operating Lamps



Cod. **ESL-50**



Cod. **ESL-60**

Cod.	ESL-50	ESL-60
Light intensity:	40000 Lux - 50000 Lux	50000 Lux - 60000 Lux
Color Temperature:	4000+/-500k	4000+/-500k
Diameter of luminous Spot:	150+/-50mm	150+/-50mm
Colour rendering index:	Ra ≥ 85	Ra ≥ 85
Power supply (VAC/DC):	220V/50Hz 110V/60Hz	220V/50Hz 110V/60Hz
Bulbs:	4 pcs	5 pcs
Rated Voltage of Bulb(V):	24V	24V
Rated Power of Bulb(W):	25W	25W
Consumption Power(W):	100W	125W
Battery:	Around 3 hour after full charge	Around 2,5 hour after full charge
Battery charge type:	Quick charge & Floating charge	Quick charge & Floating charge



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)





# ESL-50 & ESL-60

## Lampes Operatoires d'Urgence



Cod. **ESL-50**



Cod. **ESL-60**

Cod.	<b>ESL-50</b>	<b>ESL-60</b>
Intensité lumineuse:	40000 Lux - 50000 Lux	50000 Lux - 60000 Lux
Température de couleur:	4000+/-500k	4000+/-500k
Diamètre du spot lumineux:	150+/-50mm	150+/-50mm
Indice de rendu chromatique:	Ra ≥ 85	Ra ≥ 85
Alimentation (VAC/DC):	220V/50Hz 110V/60Hz	220V/50Hz 110V/60Hz
Ampoules:	4 pcs	5 pcs
Tension nominale de l'ampoule (V):	24V	24V
Puissance nominale de l'ampoule (W):	25W	25W
Consommation électrique (W):	100W	125W
Batterie:	Environ 3 heures après la charge complète	Environ 2,5 heures après la charge complète
Type de charge de la batterie:	Charge rapide & Charge d'entretien	Charge rapide & Charge d'entretien



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

