



DATA SHEET



Model: ESLD-735

Mobile Led Lamp



Description

1/ illumination	80000-120000lux
2/ color temperature	4200 \pm 300K
3/ diameter spot	250 \pm 50mm
4/ restore color index	93
5/ shadowless proportion	90%
6/ LED using life	>50000h
7/ work distance	500-2000mm
8/ up doctor head temperature	<1°C
9/ heigh adjustable	160-230cm
10/ back-up battery 12V	optional
11/ power	110-240V; 50/60Hz



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

