



EXAMINATION LED SURGICAL LIGHTS

OMEGA 200+200

Adjustable Colour Temperature

Special Features

Shadowless Light, Body is made of high quality fire resistant polymer.
Single non winged round shape lamp head for easy maintenance and better sealing against dust and fumigation.



Technical Data

Intensity	80,000 + 80,000 Lux ±10%
Size of Light Field	15 cm
Colour Temperature	5000 K
Colour Rendering Index	93 RA
LED Life	>50000 Hrs.
Number of LED	28 +28
Diameter of Light	200mm x 2
Brightness Control	Capacitive Touch Panel
Power Supply	220V/50Hz AC
Colour of LED	White
Focusing	Fixed
Battery Backup	Optional

OMEGA 200

Adjustable Colour Temperature

Special Features

Shadowless Light, Body is made of high quality fire resistant polymer.
Single non winged round shape lamp head for easy maintenance and better sealing against dust and fumigation.



Technical Data

Intensity	80,000 ±10%
Size of Light Field	15 cm
Colour Temperature	5000 K
Colour Rendering Index	93 RA
LED Life	>50000 Hrs.
Number of LED	28
Diameter of Light	200mm
Brightness Control	Capacitive Touch Panel
Power Supply	220V/50Hz AC
Colour of LED	White
Focusing	Fixed
Battery Backup	Optional

