



DATA SHEET



Model: ESK-A98/4

Electric Obstetric Bed



Description

The obstetric beds are designed for various stages of care for patients
They can be used in labor, delivery and postpartum recovery period.
The chair posture can be easily converted to provide adequate comfort for the lying-in women.
The lying-in women don't need to be moved from one bed to another thanks to the bed's perfect structure

Various actions are controlled by remote controller
Super mute big foot wheel which is easily moved and equipped with brake device.
Mattress with high density memory sponge, no seam
Length: 2150 mm
Width: 980mm
Height Adjustment: (500-900mm)
Backboard upwards: 0-65°
Trendelenburg/Rev-Trend: 7th/7th
Foot board upwards: 5°
Foot board outwards: 55°
Footboard height adjustment: 100mm

Mains voltage: AC220±10% 50HZ
Input power: 150VA
removable leg section
Weight capacity 250 kg
3 linear motors
Four 5" double side castors, central locking
ABS head board with safe lock
ABS side rails with angle indicator
Pair of hand grips
Pair of leg holders
2 IV pole forecasts
Power line
Mattress 8 cm

Standard Accessories:

A pair of leg support frame
stainless steel sewage basin
stand for transfusion
remote controller
power wire



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

